


# TwoFlex™ 2 layer Compression Bandage System

TwoFlex bandage system delivers therapeutic compression to manage venous disease and associated edema.

## APPLICATION INSTRUCTIONS

- **Apply** layer 1 first.
- **Repeat** steps 1–5 for layer 2.




**Flex the foot to 90 degrees** and apply two anchor turns starting at the base of the toes.

**Note:** Apply layer 1 with minimal tension, only the amount sufficient to cause the roll to unwind.



**Proceed up and around the ankle**, covering the Achilles tendon.




**Return over the top of the ankle and down to secure the heel.** Cover all exposed skin of the heel and ankle. Avoid creating a crease or positioning the edge of the bandage on the base of the ankle.



**Wrap in a spiral up the leg.** While wrapping, overlap the bandage by half its width.

**Note:** Apply layer 2 (short stretch cohesive bandage) at or near full stretch.



**Finish at the notch just below the knee.** Cut and discard any unused bandage. Secure with tape provided.

**Note:** Layer 2 will adhere to itself, securing the system in place. Gently press the finished wrap to the leg for best cohesion.

**Frequency of dressing change:** TwoFlex may be left in place for up to 7 days. During initial treatment, TwoFlex may need to be changed more frequently due to excessive drainage.

**Contraindications:** TwoFlex and TwoFlex XL are contraindicated for a patient with severe arterial occlusive disease with an ABI of less than 0.8. TwoFlex Lite is contraindicated for a patient with severe arterial occlusive disease with an ABI of less than 0.5.

**For more information,** see product packaging, contact your Medline Representative or call 1-800-MEDLINE (1-800-633-5463).



**Make skin health second nature**

Find more skin health insights and expertise at [medlineskinhealth.com](https://www.medlineskinhealth.com)

