



# ID NOW<sup>™</sup> COVID-19 2.0 QUICK REFERENCE INSTRUCTIONS

## Materials Required to Run a Test:



Test Base



Sample Receiver




Transfer Cartridge



Patient Swab

## Instrument:



## CLIA Complexity: WAIVED

Refer to the ID NOW COVID-19 2.0 Product Insert and User Manual for complete instructions.

## SPECIMEN COLLECTION and HANDLING

Use freshly collected specimens for optimal test performance. Inadequate specimen collection or improper sample handling may yield erroneous results.

### Nasal Swab

To collect a nasal swab sample, carefully insert the swab into the nostril exhibiting the most visible drainage, or the nostril that is most congested if drainage is not visible. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril). Rotate the swab several times against the nasal wall then slowly remove from the nostril. Using the same swab, repeat sample collection in the other nostril. Puritan® PurFlock Ultra® Flocked Swabs-Standard Tip, Copan Rayon Standard Tip Swabs, and Jiangsu Changfeng Medical Industry (JCF) Polyurethane Foam Swabs are not suitable for use. Please refer to product insert for lists of acceptable swabs for nasal swab collection.

### Nasopharyngeal Swab

To collect a nasopharyngeal swab sample, carefully insert the swab into the nostril exhibiting the most visible drainage, or the nostril that is most congested if drainage is not visible. Pass the swab directly backwards without tipping the swab head up or down. The nasal passage runs parallel to the floor, not parallel to the bridge of the nose. Using gentle rotation, insert the swab into the anterior nares parallel to the palate advancing the swab into the nasopharynx, leave in place for a few seconds, and then slowly rotate the swab as it is being withdrawn. To ensure proper collection, the swab should be passed a distance that is halfway of that from the nose to the tip of the ear. This is about half the length of the swab. **DO NOT USE FORCE** while inserting the swab. The swab should travel smoothly with minimal resistance; if resistance is encountered, withdraw the swab a little bit without taking it out of the nostril. Then elevate the back of the swab and move it forward into the nasopharynx. Puritan® Mini Rayon Tip and Puritan® PurFlock Ultra® Flocked Swabs-Mini Tip are not suitable for use in this assay. Please refer to product insert for lists of acceptable swabs for nasopharyngeal swab collection.

## Reference

Manual of Clinical Microbiology, 11th Edition, Vol. 1, ASM. (2015) pg. 279.

## Technical Support Advice Line

Further information can be obtained from your distributor, or by contacting Technical Support on:

US + 1 855 731 2288 ts.scr@abbott.com

Abbott Diagnostics Scarborough, Inc.

10 Southgate Road

Scarborough, Maine 04074 USA

www.globalpointofcare.abbott



www.globalpointofcare.abbott

## QC TESTING


External Positive and a Sterile Swab for Negative Control should be tested once with each new shipment received and once for each untrained operator.

For **QC TESTING**, select Run QC Test on the Home screen, and follow the displayed instructions. Refer to Running a QC Test in the ID NOW Instrument User Manual for further details.


1 Press Run QC Test.




2 Select COVID-19 from the menu.



3 Select the QC Test to be Run.



4 Confirm Test. Confirm the test type to match the QC sample intended for testing by pressing OK and following the on screen prompts to complete testing.



Note: The QC test is run in the same manner as a Direct Nasal Swab Patient Test. See the other side for step by step instructions.

© 2023 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Software © 2023 Axxin, used under license. All trademarks referenced are trademarks of their respective owners.

IN192001 Rev. 2 2023/09 SAP: 40003807

# ID NOW<sup>™</sup> COVID-19 2.0 QUICK REFERENCE INSTRUCTIONS

Scan to watch  
ID NOW how to  
video.



[www.globalpointofcare.abbott/en/support/product-installation-training/id-now-training-videos.html](http://www.globalpointofcare.abbott/en/support/product-installation-training/id-now-training-videos.html)


Before performing this test, refer to the ID NOW COVID-19 2.0 Product Insert and User Manual for complete Test Procedure and additional information. **Note:** To perform a quality test, refer to the **QC TESTING** section on the back.

To perform a patient test, follow directions below.

## 1 Put on a clean pair of gloves




## 2 To perform a test



- Press **Run Test**
- Select **COVID-19**
- Then follow the **ON SCREEN** instructions

## 3 Insert Test Base




## 4 Insert Sample Receiver



**!** **DO NOT REMOVE** the foil seal until instructed after 3 minute warm up.

## 5 Remove seal & mix swab




## 6 Press white Transfer Cartridge into blue Sample Receiver



- !** 1. Press **FIRMLY** the white Transfer Cartridge into blue Sample Receiver.
- 2. Listen for **CLICK(s)**.
- 3. Confirm orange indicator has risen up **FULLY** to top of cartridge.

## 7 Lift and press Transfer Cartridge into Test Base



1. Press down **FIRMLY** with **BOTH HANDS**.
2. Listen for **MULTIPLE CLICKS**.
3. Confirm orange indicator is **FULLY DOWN** before continuing.

## 8 CHECK: Make sure orange indicator is fully down



If not, **press down again**.




**Note:** Failure to fully press down may cause an **INVALID** or **FALSE** test result.

## 9 Close the lid to start test



## 10 Reading the results



Refer to Product Insert instructions for retest procedure in the case of an Invalid result.


## 11 Optional workflow to run an Influenza A & B test

### FOR SEQUENTIAL WORKFLOW (COMBO):

If Sequential Workflow is desired.

**DO NOT REMOVE** Sample Receiver!

- Press **Actions**
- Select **COMBO: COVID-19 + Flu A & B**





**Note:** See Sequential Workflow Instruction in Product Insert IN19200 for details.

## 12 Discard used test pieces



Refer to the Product Insert for instructions on safe handling and disposal of samples and test components.

<p>Abbott ID NOW COVID-19 2.0</p> <p>QRI, EN</p> <p><b>Size:</b> 8.5 in x 14 in</p>	<p style="text-align: center;"><b>Printed Colors</b></p>  <p style="text-align: center;">CMYK</p> <hr/> <p style="text-align: center;"><b>Incoming Inspection Colors (For Reference Only)</b></p> <p style="text-align: center;"><small>Colors below are not used for printing</small></p>  <p style="text-align: right;"><small>PMS 2995 U Primary Blue</small></p> <p style="text-align: right;"><small>PMS 303 U Dark Blue</small></p> <p style="text-align: right;"><small>PMS 224 U Magenta-Pink</small></p> <p style="text-align: right;"><small>PMS 185 U Red</small></p>	<p>PN: IN192001 Rev: 2 SAP: 40003807</p> <p><b>Date of Last Revision:</b> 2.8 2023/09/14</p>
-------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------