

# MCKESSON

## Raised Toilet Seat removable padded arms plastic frame

Provides a total weight bearing capacity of 300 lbs (136 kg).

For use with MFR#'s 146-RTL12027RA.

### INSTALLATION INSTRUCTIONS

1. Remove product from carton and READ ALL INSTRUCTIONS CAREFULLY before attempting to install product.
2. To install arms, depress push pins on arm and insert into mounting tubes. There will be an audible click. Ensure the push pins fully protrude through the holes in the mounting tubes. Make sure arms are locked in place by pulling them. To remove arms, reverse steps.
3. Lift toilet lid and seat assembly and place rear wings of toilet seat around rim of toilet bowl. (Some toilet bowls may require that seat assembly be removed to allow proper fit.) In addition the front adjustment bracket must clear the rim in order for the raised toilet seat to fit securely.
4. With raised toilet seat resting correctly on rim, straighten and center product over toilet bowl as necessary. Turn adjusting knob clockwise to tighten product securely in place. Hand tighten only - use no tools, and be careful not to over-tighten, as threads on adjustment bracket may strip.
5. Check for secure fit routinely before using product. User's weight should be centered over toilet seat, not too far forward or too far to one side, as unit may tip.
6. If you have a style of commode where this product will not fit properly, or cannot be tightened securely as noted above, it is recommended that you discontinue usage.

### REMOVAL INSTRUCTIONS

1. Rotate knob counter clockwise until seat is loose.
2. Lift front of raised toilet seat and pull forward to move.

### CLEANING INSTRUCTIONS

#### CAUTION

1. Always be sure that the seat is correctly and securely locked in place before using.
2. Be sure that the adjustment knob in front of the seat is tight at all times, thus securing the seat to the toilet bowl.
3. This product is designed for use with the most common types of commodes. It will not fit all commode types; be sure that the unit is correctly installed, fits properly and is stable before use.
4. On models with arm rests: The arms are for the purpose of providing assistance when getting on and off the seat. DO NOT try to use the arms to support full body weight during transfer, as this may result in failure of the unit.

Technical support, replacement parts and warranty information? Call 1-866-694-5085. General product questions? Call 1-800-777-4908

#### Satisfaction Guaranteed

For complete details, please visit  
[mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

Distributed by McKesson Medical-Surgical Inc.  
Richmond, VA 23233  
PVN B0518  
Made in China