

McKESSON

Low Suds Liquid Instrument Detergent

FOR EFFECTIVE CLEANING OF MEDICAL DEVICES AND SURGICAL INSTRUMENTS

For use with MFR #'s 53-28551

McKesson Low Suds Liquid Instrument Detergent is designed to effectively clean medical devices and surgical instruments. Its low foaming, neutral pH formula is appropriate for use in automated washers as well as manual cleaning systems. **McKesson Low Suds Liquid Instrument Detergent** is non-corrosive, contains no chlorines, phosphates, or hydroxides, and will not harm surgical instruments.

DIRECTIONS FOR USE:

1. Always wear gloves and eye protection.
2. Pour product into basin or sink at a dilution of 1/8 ounce to 2 ounces (3.7 mL to 59 mL) of product per one gallon of water.
3. Dip instruments into solution, clean and remove.
4. Rinse instruments in water.
5. Wrap instruments and proceed with further processing.
6. When using automated washers, always consult equipment manufacturer guidelines for recommended dilutions and other pertinent information. Make certain machine goes through both wash and rinse cycles.

STORAGE

Store in a dry place no lower in temperature than 50°F (10°C) or higher than 120°F (49°C).

CAUTION

KEEP OUT OF REACH OF CHILDREN

Eyes: Do not get into eyes. Avoid contact with eyes by wearing eye protection. In case of eye contact, flush immediately with plenty of water. Get medical attention.

Skin: Avoid skin contact. Prolonged or repeated contact may cause irritation or skin sensitization. In case of skin contact, rinse with water. Get medical attention if irritation occurs.

Ingestion: May be harmful if swallowed. If swallowed, drink large quantities of water or milk. **DO NOT INDUCE VOMITING.** Do not give anything by mouth to an unconscious person. Get medical attention.

Questions? Call 1-800-777-4908

■ Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Distributed By McKesson Medical-Surgical Inc.
Richmond, VA 23228
Rev. 00 01/23
Made in USA