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
## Instructions for use

### General description of the device

The TwoPress 2 Lite\* 2-component compression system is a combination of two non-woven bandages. One of the bandages is made of polypropylene and elastane and is brown in color, the other one is made of polyester and elastane and is white in color. The brown one has a cohesive layer on both sides, the white one has a cohesive layer on one side. The non-active and non-sterile devices are available as a kit. Both bandages are 10 cm in width. TwoPress 2 Lite contains indicator prints on both layers.

### Intended purpose

The 2-component compression system is a single-use, non-sterile medical device used by healthcare professionals in clinical and home environment. It can be used for compression therapy which requires a reduced level of compression (ABPI 0.6-0.8), or for patients unable to tolerate full compression, or as part of a step-up approach, for the treatment of:

- mixed etiology (arterial-venous) leg ulcers
- venous leg ulcers
- edema
- thrombosis/thrombophlebitis

The 2-component compression system can be used for long-term compression therapy on intact skin in humans (above 6 months of age). It can be used in combination with primary and secondary dressings.

### Composition

Short-stretch padding bandage: 92% polyester + 8% elastane, with a cohesive coating on one side.

Compression bandage: 85% polypropylene + 15% elastane, with a cohesive coating on both sides.

### Indication

The 2-component compression system can be used for compression therapy which requires a reduced level of compression (ankle-brachial pressure index [ABPI] 0.6-0.8), or for patients unable to tolerate full compression, or as part of a step up approach, for the treatment of:

- mixed etiology (arterial-venous) leg ulcers
- venous leg ulcers
- edema (e.g. venous edema, lymph edema, post-traumatic edema)
- thrombosis/thrombophlebitis (e.g. thrombophlebitis as well as condition after healed thrombophlebitis, deep leg and arm vein thrombosis, condition after deep vein thrombosis, post-thrombotic syndrome, thromboprophylaxis)

### Intended users

TwoPress 2 Lite is designed to be used by healthcare professionals who have received appropriate training.

### Contraindications

The 2-component compression system must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD) with critical ischemia
- decompensated heart failure (NYHA III + IV)
- phlegmasia coerulea dolens
- ankle-brachial pressure index (ABPI) < 0.6
- septic phlebitis
- florid erysipelas

### Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral neuropathy)
- pronounced weeping dermatoses
- incompatibility with compression bandages material
- primary chronic polyarthritis

The application technique, e.g. with regard to stretch and/or padding, may have to be adapted to the patient's anatomical and pathophysiological characteristics, especially in case of fragile or elderly skin. Wounds have to be covered adequately.

If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. The right measurement method must be considered for the ABPI. During the use of the 2-component compression system, the ABPI must be controlled regularly. The therapy must be adapted accordingly.

The devices should not be used on infants below the age of 6 months. When applying the 2-component compression system, on children (from 6 months on), the pressure should be reduced appropriately.

If at any time during therapy the patient develops severe pain, numbness, tingling or noticeable discoloration of the toes, a healthcare professional must be contacted immediately, and the 2-component compression system must be removed.

### Application

Please use this system only as a set and in the correct order. When applying the two bandages, make sure that the patient's foot is at an angle of 90° relative to the lower leg (neutral 0 position) ❶. Examine the leg before application to identify bony prominences and potential pressure points and pad them with the short-stretch padding bandage. Start applying the white short-stretch padding bandage (1) with the non-cohesive side against the skin proximal to the base of the toes. After two turns of fixation with adequate pressure (depending on the individual foot position (pronation or supination)), the placement of the bandage should be chosen medially or laterally). Wrap the bandage across the back of the foot and enclose the heel in figures of 8. Wrap the foot by pressing adequately while considering anatomical / non-cylinder shape (according to Laplace's law) which, for instance, creates a pressure increase laterally (small radius) and a pressure decrease on the forefoot (large radius). Ignore the printed indicator so far. ❷❸. On the lower leg, wrap the white short-stretch padding bandage around it using a spiral technique with an overlap of 50% and following the indicator for optimal stretch ❹. To achieve the optimal pressure, stretch the bandage until the indicator that is printed on it has reached the length of a regular hexagon with all sides having the same length. Make sure that the bandage is not overstretched. Pull the bandage forward before every half-turn and place the bandage with minimal pull over the mid-line of the printed hexagon below. After each turn, place the green pressure indicator on the printed hexagons so that both can match. The stretch may have to be adapted to the patient's anatomical and pathophysiological characteristics. Tear off any excess bandage material beneath the tuberositas tibiae and about 2 cm below the fibular head ❺. Additional fixation is not necessary; tape can be used if required.

As with the white short-stretch padding bandage, start applying the brown cohesive compression bandage (2) proximal to the base of the toes, so that the white short-stretch padding can still be seen here. After two turns for fixation, wrap the bandage across the back of the foot and enclose the heel in figures of 8. Wrap the foot by pressing adequately while considering anatomical / non-cylinder shape (according to Laplace's law) which, for instance, creates a pressure increase laterally (small radius) and a pressure decrease on the forefoot (large radius). Ignore the printed indicator so far. ❶. Wrap the brown cohesive compression bandage using a spiral technique around the lower leg with an overlap of 50% ❷. To achieve the optimal pressure, stretch the bandage until the indicator that is printed on it has reached the length of a regular hexagon with all sides having the same length. Make sure that the bandage is not overstretched. Pull the bandage forward before every half-turn and place the bandage with minimal pull over the mid-line of the printed hexagon below. After each turn, place the green pressure indicator on the printed hexagons so that both can match ❸. The overlap must be 50%. For better orientation, use the center line of the hexagon ❹. Check the elongation with green pressure indicator to check also the correct stretch of the second layer. If the indicator is correctly stretched, you will see the printed indicator at the edge of the template. The stretch may have to be adjusted to the patient's anatomical and pathophysiological characteristics. Tear off any excess bandage material beneath the tuberositas tibiae and about 2 cm below the fibular head, so that the white short-stretch padding is still visible. For optimal cohesion, press the compression bandage gently and evenly with both hands over the entire lower leg ❺❻. Additional fixation is not necessary; tape can be used if required.

### Further remarks

TwoPress 2 Lite is indicated for use by medical professionals who have received appropriate training. Before applying TwoPress 2 Lite, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. foot, anterior crest of tibia, achilles tendon, shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high. TwoPress 2 Lite is intended for single use only. Both bandages can be torn by hand. Both bandages are made without natural rubber latex. The applied bandage can be removed either by unwrapping, or by carefully cutting it with bandage scissors. The bandage can be worn day and night for up to 7 days. The specific application period must be determined individually based on the particular patient's wound/skin/edema situation. The application technique must be followed in order to achieve the therapeutically recommended compression pressure. If the bandage becomes undone or slips, the bandage should be removed, and a new system applied in consultation with your medical staff.

If at any time during therapy with TwoPress 2 Lite the patient develops severe pain, numbness, tingling or noticeable discoloration of the toes, a medical professional must be contacted immediately, and the bandage must be removed. When applying TwoPress 2 Lite, make sure that the tension is even as the hexagon indicator reaches the shape of a regular hexagon ❹ and the bandages are straight in order to achieve uniform pressure distribution and avoid constrictions. During therapy with TwoPress 2 Lite, the patient should be encouraged to move. An important element of compression therapy is the active muscle pump, which acts as counterpart to the working pressure of the compression bandage.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices), if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

### Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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\* © MEDICE-Iserlohn