





Extension Feeding Set Attachment Instructions and Troubleshooting Tips

Attachment to Gastric Feeding Device

Always be sure to thoroughly wash hands before handling feeding devices. Before attaching the feeding set, it is important to verify the correct placement of the G-Button or other gastric feeding device inside the stomach and its functionality as defined in the instructions for use of the G-Button or other gastric feeding device.



INSTRUCTIONS	TROUBLESHOOTING TIPS	ILLUSTRATIONS
<p>1. Remove the cap from the G-Button or other gastric feeding device to reveal the anti-reflux valve. (Figure 2).</p>	<p>Steady the G-button by holding the feeding port with thumb and index finger with one hand (Figure 3A or 3B).</p>	
<p>2. Connect the extension feeding set's secure lock connector (feeder) to the G-button or other feeding device. (Figure 2).</p>	<p>Make sure the feeder and connecting tooth are directly aligned to feeding port opening (Figure 3A or 3B).</p>	
<p>3. Align the connecting tooth to the groove in the connecting ring. Once inside, press gently on the feeder and turn half a turn, clockwise, to secure the feeder in place.</p>	<p>Gently press feeder with one hand while holding the feeding port with the other.</p>	
<p>CAUTION: If the feeder is not properly connected and secure, leakage may occur.</p>	<p>Make sure the connecting tooth is not visible and the feeding set is locked in place to avoid leakage.</p>	
<p>4. Confirm that the device is connected to an enteral port and not an IV set.</p>	<p>Check that the connection between the feeding set and the G-button has a secure fit like the illustration.</p>	
<p>5. Flush the extension set with warm water (10mL for adults or 3-5mL for children).</p>		
<p>6. Check for ease of flushing.</p>	<p>Water for flushing should flow freely with only gentle thumb pressure applied to the syringe.</p>	
<p>After decompression or enteral feeding, flush the feeding set with 5-10mL of warm water.</p>		
<p>Removing Feeding Set</p>		
<p>Gently turn the feeder counter-clockwise until the connecting tooth can be released through the groove.</p>	<p>Once the feeding set is properly removed, be sure to follow instructions for proper washing and reuse.</p>	

Extension Feeding Set Attachment Instructions and Troubleshooting Tips

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.