Fully Automatic Arm Cuff Blood Pressure Monitor

Electronic Sphygmomanometer

Operation Guide

REF MDS3001 REF MDS3001U

REF MDS3001PLUS

CONTENTS

Important Information	3	Displaying Stored Results	22
Contents And Display Indicators	4	Deleting Measurements	
Intended Use	5	From The Memory	24
Contraindication	5	Assessing High Blood Pressure For Adult	24
Product Description	6	Medline Blood Pressure Cuff Table	26
Specifications	8	Troubleshooting	27
Notice	10	Maintenance	30
Setup And Operating		Explanation Of Symbols On Unit	32
Procedures	15	Warranty Information	33
Battery Loading and AC Adapter Loading	15	Service Center	33
Clock And Date Adjustment	16	Electromagnetic	
Technical Alarm Description	17	Compatibility Information	34
Connecting The Cuff To The Monitor	17		
Applying The Cuff	18		
Body Posture During Measurement	19		
Taking Your Blood Pressure Reading	20		

IMPORTANT INFORMATION

NORMAL BLOOD PRESSURE FLUCTUATION

All physical activity, excitement, stress, eating, drinking, smoking, body posture and many other activities or factors (including taking a blood pressure measurement) will influence blood pressure value. Because of this, it is unusual to obtain identical multiple blood pressure readings.

Blood pressure fluctuates continually day and night. The highest value usually appears in the daytime and the lowest one usually at midnight. Typically, the value begins to increase at around 3:00 AM, and reaches the highest level in the daytime while most people are awake and active.


Considering the above information, it is recommended that you measure your blood pressure at approximately the same time each day.

Please always relax a minimum of 1 to 15 minutes between measurements to allow the blood circulation in your arm to recover. It is unusual to obtain identical blood pressure readings each time.

Too frequent measurements may cause injury due to blood flow interference.

CONTENTS AND DISPLAY INDICATORS

· 1 Blood Pressure Monitor · 1 Operation Guide · 1 Arm Cuff · 1 Soft Storage Case

Picture 1

INTENDED USE

Fully Automatic Electronic Blood Pressure Monitor is for use at home and is a non-invasive blood pressure measurement system intended to measure the diastolic and systolic blood pressures and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm.

CONTRAINDICATION

 It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

PRODUCT DESCRIPTION

Based on Oscillometric methodology and silicon integrated pressure sensor, blood pressure and pulse rate can be measured automatically and non-invasively. The LCD display will show blood pressure and pulse rate. The most recent 90 measurements can be stored in the memory with date and time stamp. The monitor can also show the average reading of the last three measurements. The Electronic Sphygmomanometer corresponds to the below standards: IEC 60601-1Edition 3.1 2012-08/EN 60601-1: 2006/A1: 2013 (Medical electrical equipment – Part 1: General requirements for basic safety and essential performance), IEC60601-1-2: 2014/EN 60601-1-2: 2015 (Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC80601-2-30: 2018/EN IEC80601-2-30: 2019 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) EN 1060-1: 1995 + A2: 2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3:

1997 + A2: 2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems);
ISO81060-2: 2013 (Non-Invasive Sphygmomanometers - Part 2: Clinical Validation Of Automated Measurement Type).

SPECIFICATIONS

1. Product name: Blood Pressure Monitor
2. Model: MDS3001 / MDS3001PLUS / MDS3001U
3. Classification: Internally powered, Class II, Type BF applied part, IP20, No AP or APG, Continuous operation
4. Machine size: 165mmx96mmx65mm (6 1/2"x 3 25/32"x 2 9/16")
5. Cuff circumference:
Standard Adult Cuff 22cm - 30cm (8 21/32" - 11 13/16"),
Large Adult Cuff 30cm - 42cm (11 13/16" - 16 17/32") ,
Extra Large Adult Cuff 42cm - 48cm (16 17/32"- 18 29/32")
Universal Cuff 22cm - 42 cm (8 21/32" - 18 29/32")
6. Weight: approx. 300g (10 19/32 oz.) (excluding batteries and cuff)
7. Measuring method: oscillometric method, automatic air inflation and measurement
8. Memory volume: 90 times with time and date stamp
9. Power source: Batteries: 4 x1.5V = SIZE AA
AC adapter:
 - Input: 100-240V ~ 50/60Hz, 60mA
 - Output: DC 6.0V = 600mA
10. Measurement range: Cuff pressure: 0-300mmHg

Systolic: 60-260mmHg
Diastolic: 40-199mmHg
Pulse rate: 40-180 beats/minute

11. Accuracy: BP Measurement: This monitor is validated to meet the accuracy requirements of ± 5 mmHg as indicated in the clinical standard(ISO81060-2)
Cuff Pressure/Sensor: ± 3 mmHg
Pulse Rate: 40-180 beats/minute
Less than 60: ± 3 beats/minute
More than 60: $\pm 5\%$
12. Environmental temperature for operation: $10^{\circ}\text{C}\sim 40^{\circ}\text{C}$ ($50^{\circ}\text{F}\sim 104^{\circ}\text{F}$)
13. Environmental humidity for operation: $\leq 85\%$ RH
14. Environmental temperature for storage and transport: $-20^{\circ}\text{C}\sim 55^{\circ}\text{C}$ ($-4^{\circ}\text{F}\sim 131^{\circ}\text{F}$)
15. Environmental humidity for storage and transport: $\leq 90\%$ RH
16. Environmental pressure: 80KPa-105KPa
17. Battery life: Approx 540 times.
18. List of all components belonging to the pressure measuring system, including accessories: Pump, Valve, LCD, Cuff, Sensor

Note: These specifications are subject to change without notice.

NOTICE

1. Stay still, remain calm and rest for 5 minutes before blood pressure measurement.
2. The cuff should be placed at the same level as your heart.
3. During measurement, neither speak nor move your body and arm.
4. Measure on the same arm each time.
5. Please always relax a minimum of 1 to 15 minutes between measurements to allow the blood circulation in your arm to recover. Prolonged over-inflation (cuff pressure exceeding 300 mmHg or maintained above 15 mmHg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.
6. Remove the batteries if the monitor will not be used for a month or more to avoid damage caused by battery leakage.
7. ⚠️ This Electronic Sphygmomanometer is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
8. Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard Institute, electronic or automated sphygmomanometers.

9. Information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference please see part ELECTROMAGNETIC COMPATIBILITY INFORMATION. It is suggested that the blood pressure monitor be kept at least 30 cm away from other wireless devices, such as WLAN unit, microwave oven, etc. It can't be used near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high
10. If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of "♥" will be displayed. Under this condition, the Electronic Sphygmomanometer can continue functioning, but the results may not be accurate; it's suggested that you consult with your physician for accurate assessment. There are 2 conditions under which the signal of IHB will be displayed:
- a) The coefficient of variation (CV) of pulse period >25%.
 - b) The difference of adjacent pulse period $\geq 0.14s$, and the number of such pulse takes more than 53 percent of the total number of pulse.
11. Read all the information in the operation guide and any other literature in

- the box before operating the unit.
12. Please only use the cuff that was supplied by the manufacturer; otherwise a biocompatibility hazard can occur, resulting in a measurement error.
 13. ⚠ The monitor might not meet its performance specifications or cause a safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.
 14. ⚠ Please do not share the cuff with infected people to avoid cross-infection.
 15. Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
 16. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the

equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

17. Medical AC adapter which input is 100-240V ~ 50/60Hz, 60mA and output is DC 6.0V, 600mA and complied with IEC 60601-1/EN60601-1/UL 60601-1 and IEC 60601-1-2/EN 60601-1-2/UL 60601-1-2 is suitable for this monitor. Please note that the monitor jack size: hole \varnothing 5.5mm, center pin \varnothing 2.0mm. Please pay attention to polarity.

18. Consult with your physician for any of the instances below.

- This device does not measure circulation or blood flow.
- The application of the cuff is over a wound or inflammation.
- The application limb shares intravascular access or therapy, or an arteriovenous (A-V) shunt.
- The application arm is on the same side as a mastectomy or lymph node clearance.
- The application limb shares monitoring equipment.

19. The patient is an intended operator.
20. Attention that changes or modification not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
21. Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the batteries and the unit out of the reach of children and disabled persons.
22. If you are allergic to plastic/rubber, please don't use this device.
23. The device is not intended for use on neonates, children or pregnant women. (Clinical testing has not been conducted on neonates, children or pregnant women.)
24. a) Motion, trembling, shivering may affect the measurement reading.
b) Consult your physician before using the device for any of the following conditions: common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pre-eclampsia, renal diseases.
25. Do not use this unit in a moving vehicle, this may result in erroneous measurement.

SETUP AND OPERATING PROCEDURES

1. BATTERY LOADING AND AC ADAPTER LOADING

- a. Open battery cover at the back of the monitor.
- b. Load four “AA” size batteries. Please pay attention to polarity.
- c. Close the battery cover.
 - ⚠ When LCD shows battery symbol , replace all batteries with new ones.
 - ⚠ Rechargeable batteries are not suitable for this monitor.
 - ⚠ Remove the batteries if the monitor will not be used for a month or more to avoid significant damage caused by battery leakage.
 - ⚠ Avoid getting the battery fluid in your eyes. If it should get in your eyes, immediately rinse with plenty of clean water and contact a physician.
 - ~~⚠~~ The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.
- d. If you use the AC adapter, please make sure that the monitor is turned off and the batteries have been removed. Put the connector plug of the adapter into the monitor socket. Then plug the adapter to AC electrical outlet. To disconnect the AC adapter, remove the AC adapter from the electrical outlet, and remove the AC adapter plug from the monitor socket.

Please note: AC adapter is not included.

Do not plug or unplug the power cord into the electrical outlet with wet hands.

Do not overload power outlets. Plug the device into the appropriate voltage outlet. If the AC adapter isn't working, please change the adapter.

Do not pull out the adapter when you are using the monitor.

Do not use any other type of AC adapter as it may harm the monitor.

Do not to position the monitor so that it is difficult to operate the disconnection device.

Picture 2

2. CLOCK AND DATE ADJUSTMENT

- a. Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date. See picture 2.
- b. While the monitor is in Clock Mode, press both the "START" and "MEM" button down, then release the buttons. The month will blink at first. Press the button "START" repeatedly, the day, hour and minute will blink in turn. While the number is blinking, press the button "MEM" to increase the number. Keep on pressing the button "MEM" to increase the number fast.
- c. You can turn off the monitor by pressing "START" button when the minute

- is blinking, then the time and date is confirmed.
- d. The monitor will turn off automatically after 1 minute of no operation.
 - e. Once you change the batteries, you should readjust the time and date.

3. TECHNICAL ALARM DESCRIPTION


The monitor will show 'HI' or 'Lo' as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICATIONS. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and does not require a reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

4. CONNECTING THE CUFF TO THE MONITOR

Insert the cuff tubing connector into the socket in the left side of the monitor. Make certain that the connector is completely inserted to avoid air


Picture 3

leakage during blood pressure measurements.

- ⚠️ Avoid compression or restriction of the connection tubing during measurement, which may cause inflation error, or harmful injury due to continuous cuff pressure.

5. APPLYING THE CUFF

- Pulling the cuff end through the metal loop (the cuff is packaged like this already), turn it outward (away from your body) and tighten it and close the velcro fastener.
- Place the cuff around a bare arm 1-2 cm above the elbow joint.
- While seated, extend your arm with the palm facing up in front of you on a flat surface such as a desk or table. Position the air tube in the middle of your arm in line with your middle finger.
- The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger between your arm and the cuff.

Picture 4

Note:

- Please refer to the cuff circumference range in “SPECIFICATIONS” to make**


sure that the appropriate cuff is used.

- **Measure on the same arm each time.**
- **Do not move your arm, body, or the monitor and do not move the tube during measurement.**
- **Stay still, remain calm and rest for 5 minutes before blood pressure measurement.**
- **Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clean it by hand in a mild detergent, then rinse it thoroughly in cold water.**
- **Never dry the cuff in clothes dryer or iron it. Clean the cuff after 200 uses.**
- **Do not place the cuff around your arm if the arm has any inflammation, acute diseases, infections or skin wounds.**

6. BODY POSTURE DURING MEASUREMENT

Sitting Measurement


- a. Sit in a chair with your feet on the floor, do not cross your legs.
- b. Extend your arm with the palm facing up in front of you on a flat surface such as desk or table.
- c. The cuff should be at the same level as your heart.

Picture 5

Lying Down Measurement

- Lie on your back.
- Place your arm straight along your side with your palm facing down.
- The cuff should be placed at the same level as your heart.

Picture 6

7. TAKING YOUR BLOOD PRESSURE READING

Picture 7


Picture 7.1


Picture 7.2

Picture 7.3

Picture 7.4

Picture 7.5

- After applying the cuff with your body in a comfortable position, press the "START" button. A beep is heard and all display characters are shown for self-test. See picture 7. Please contact the service center if segment is missing.
- The current memory bank (U1, U2 or U3) should be blinking. See picture 7.1.


Press “MEM” button to change over to other bank. Confirm your selection by pressing “START” button. The current bank can also be confirmed automatically after 5 seconds with no operation.

- c. After selecting the memory bank, the last result will be displayed with date and time. See picture 7.2. LCD will show picture 7.3 with current date and time if no memory stored in the selected bank.
- d. The monitor will start to seek zero pressure. See picture 7.4.
- e. The monitor inflates the cuff until sufficient pressure has built up for a measurement. Then the monitor slowly releases air from the cuff and conducts the measurement. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen. Irregular heartbeat symbol (if any) will blink. See picture 7.5. The result will be automatically stored in the current memory bank.
- f. After measurement, the monitor will turn off automatically after 1 minute of no operation. Alternatively, you can press the “START” button to turn off the monitor manually.
- g. During measurement, you can press the “START” button to turn off the monitor manually.

Note: Please consult a health care professional for interpretation of pressure measurements.

8. DISPLAYING STORED RESULTS

- a. While the monitor is in Clock Mode, press “MEM” button to display the stored results. The memory bank will blink and the amount of results in this bank will be displayed. See picture 8. Press “START” button to change over to other bank. Confirm your selection by pressing “MEM” button. The current bank can also be confirmed automatically after 5 seconds with no operation.
- b. After selecting the memory bank, the LCD will display the average value of the last three results in this bank, See picture 8.1. If no result stored, LCD will show picture 8.2.
- c. Two seconds later, the most recent result will be displayed with date and time. See picture 8.3. Irregular heartbeat symbol (if any) and the indicator of blood pressure classification will blink at the same time. Press “MEM” button again to review the next result. In this way, the results measured previously in the current bank can be displayed.
- d. When displaying the stored results, the monitor will turn off automatically after 1 minute of no operation. You can also press the button “START” to turn off the monitor manually.

Picture 8


Picture 8.1

Picture 8.2

Picture 8.3


Picture 8.4

9. DELETING MEASUREMENTS FROM THE MEMORY

When any result is displayed, continuously press the “MEM” button for three seconds; all results in the current bank will be deleted after three “beeps”. LCD will show picture 9, then come to the mode without result. See picture 9.1.

Picture 9


Picture 9.1

10. ASSESSING HIGH BLOOD PRESSURE FOR ADULTS

The following guideline for assessing high blood pressure (without regard to age or gender) has been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

Classification of blood pressure for adults

BLOOD PRESSURE CLASSIFICATION	SBP mmHg	DBP mmHg	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High-Normal	130-139	85-89	GREEN
Grade 1 Hypertension	140-159	90-99	YELLOW
Grade 2 Hypertension	160-179	100-109	ORANGE
Grade 3 Hypertension	≥ 180	≥ 110	RED

WHO/ISH Definitions and Classification of Blood Pressure Levels

Medline Blood Pressure Cuff Table

If you need a different sized cuff, please consult the table below for information about the different cuff sizes Medline offers for this monitor. Contact Medline Industries, LP at 1-800-MEDLINE or visit www.medline.com if you need to purchase a new cuff. Provide the appropriate item number from the table below.

MEDLINE ITEM NUMBER	SIZE	MEASUREMENT RANGE (cm)
MDS9971	Adult	22-30 cm
MDS9972	Large Adult	30-42 cm
MDS9973	Extra Large Adult	42-48 cm
MDS9974	Universal	22-42 cm

BLOOD PRESSURE MONITOR ITEM NUMBER	INCLUDED ARM CUFF ITEM NUMBER	ARM CUFF SIZE	MEASUREMENT RANGE (cm)
MDS3001	MDS9971	Adult	22-30 cm
MDS3001U	MDS9974	Universal	22-42 cm
MDS3001PLUS	MDS9971 & MDS9972	Adult & Large Adult	22-30 cm & 30-42 cm

TROUBLESHOOTING 1

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again
	Body posture was not correct during testing	Review the "BODY POSTURE DURING MEASUREMENT" sections of the instructions and re-test
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heartbeat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer

TROUBLESHOOTING 2

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows low battery symbol 	Low Battery	Change the batteries
LCD shows "Er 0"	Pressure system is unstable before measurement	Don't move and try again
LCD shows "Er 1"	Fail to detect diastolic pressure	
LCD shows "Er 2"	Fail to detect systolic pressure	
LCD shows "Er 3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again
LCD shows "Er 4"	Pneumatic system leakage or cuff is too loose during inflation	

TROUBLESHOOTING 2.2

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows "Er 5"	Cuff pressure above 300mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory
LCD shows "Er 6"	More than 3 minutes with cuff pressure above 15 mmHg	
LCD shows "Er 7"	EEPROM accessing error	
LCD shows "Er 8"	Device parameter checking error	
LCD shows "Er 9"	MCU self-verify error	
LCD shows "Er A"	Pressure sensor parameter error	
No response when you press button or load battery	Incorrect operation or strong electromagnetic interference	Take out batteries for five minutes, and then reinstall all batteries.

MAINTENANCE

1. Do not drop this monitor or subject it to strong impact.
2. Avoid high temperature and sunlight. Do not immerse the monitor in water as this will result in damage to the monitor.
3. If this monitor is stored near freezing, allow it to acclimate to room temperature before use.
4. Do not attempt to disassemble this monitor.
5. If you do not use the monitor for a long time, please remove the batteries.
6. It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.
7. If the monitor becomes dirty, please clean it with a soft dry cloth. Do not use any abrasive or volatile cleaners.
8. No component in the monitor can be maintained by the user. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment that are designated as reparable can be supplied.
9. The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years.

10. The cuff integrity is maintained after 1,000 open-close cycles of the closure.
11. The cuff will maintain performance characteristics for a minimum of 10,000 inflations.
12. Not servicing/maintenance while the monitor is in use.
13. The monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20 °C.
14. The monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20 °C.

EXPLANATION OF SYMBOLS ON UNIT

Symbol for “THE OPERATING GUIDE MUST BE READ”(The sign background color: blue. The sign graphical symbol: white)

Symbol for “CAUTION”

Symbol for “Cuff Type: BF Applied Part”

Symbol for “ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local Authority or retailer for recycling advice”.

Symbol for “KEEP DRY”

Symbol for “POLARITY OF D.C. POWER CONNECTOR”

IP20

The first characteristic numeral symbol for “Degrees of protection against access to hazardous parts and against solid foreign objects”. The second characteristic numeral symbol for “Degrees of protection against ingress of water”. The bigger the number is, the higher the level of protection is.

Symbol for “SERIAL NUMBER”

MR Unsafe

WARRANTY INFORMATION

The Medline automatic electronic blood pressure monitor MDS3001/MDS3001PLUS/MDS3001U is under warranty for 2 years from the date of purchase. This warranty includes the instrument and the cuff. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument. The warranty is only valid upon presentation of proof of purchase.

SERVICE CENTER

Manufactured for: Medline Industries, LP,
Three Lakes Drive, Northfield, IL 60093 USA.
Tel: 1-800-MEDLINE V10 RJ24AXA

ELECTROMAGNETIC COMPATIBILITY INFORMATION

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration - Emission

The [MDS3001/MDS3001PLUS/MDS3001U] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3001/MDS3001PLUS/MDS3001U] should assure that it is used in such an environment.

Phenomenon	Compliance	Electromagnetic environment - guidance
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

Guidance and manufacturer's declaration - Enclosure Port

The [MDS3001/MDS3001PLUS/MDS3001U] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3001/MDS3001PLUS/MDS3001U] should assure that it is used in such an environment.

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to Proximity fields from RF wireless communications equipment
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Guidance and manufacturer's declaration - Proximity fields from RF wireless communications equipment

The [MDS3001/MDS3001PLUS/MDS3001U] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3001/MDS3001PLUS/MDS3001U] should assure that it is used in such an environment.

Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ± 5 kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		

810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Guidance and manufacturer's declaration - Input a.c. power Port

The [MDS3001/MDS3001PLUS/MDS3001U] is intended for use in the electromagnetic environment specified below. The customer or the user of the [MDS3001/MDS3001PLUS/MDS3001U] should assure that it is used in such an environment.

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrical fast transients/ burst	IEC 61000-4-4	±2 kV 100kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV
Surges Line-to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80%AM at 1kHz

Voltage dips	IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0% UT; 250/300 cycles

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