



HANDLED

WITH CARE.

Medline Heel Incision Devices
Safe, Easy-to-Use Heel Lancets

DESIGNED WITH THE PRECISION YOU DEMAND.



Medline Heel Incision Devices

Our easy-to-use trigger activates the blade quickly, safely and precisely making sure that each incision is made at the depth and width required for consistent blood flow.

- » High-quality stainless steel blade
- » Small, precise incisions reduce bruising
- » Color-coded for easy size identification
- » Disposable, one-time use

For optimal results, use the Medline Infant Heel Warmer to help increase blood flow to the heel.

Increase Blood Flow

To ensure consistent blood flow, prep the baby's heel with an Infant Heel Warmer. Infant Heel Warmers help minimize potential bruising, infection or cartilage damage caused by frequent heel sticks from blood tests.



- » Specially designed gel pack conforms gently to infant's heel
- » Lift-and-stick tape for easy application and removal
- » Activates instantly

Item No.	Description	Pkg.
MDS138007	Infant Heel Warmer, Gel Pack	100/cs

Item No.	Description	Pkg.
MPHMPREEMIE2	Micro-Preemie Heel Lancet, Light Blue, 0.65 mm Depth, 1.4 mm Length	200/cs
MPHMPREEMIE2Z	Micro-Preemie Heel Lancet, Light Blue, 0.65 mm Depth, 1.4 mm Length	50/bx
MHPREEMIE2	Preemie Heel Lancet, White, 0.85 mm Depth, 1.75 mm Length	200/cs
MHPREEMIE2Z	Preemie Heel Lancet, White, 0.85 mm Depth, 1.75 mm Length	50/bx
MPHNEWBORN2	Newborn Heel Lancet, Pink, 1.0 mm Depth, 2.5 mm Length	200/cs
MPHNEWBORN2Z	Newborn Heel Lancet, Pink, 1.0 mm Depth, 2.5 mm Length	50/bx
MPHTODDLER2	Toddler Heel Lancet, Purple, 2.0 mm Depth, 3.0 mm Length	200/cs
MPHTODDLER2Z	Toddler Heel Lancet, Purple, 2.0 mm Depth, 3.0 mm Length	50/bx

The defined contact area allows for accurate blade placement.



Our integrated safety lock protects against accidental finger sticks as well as eliminates any loose pieces of plastic.



Contact your Medline Representative or call 1-800 MEDLINE to order yours today.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)