

March 2023

McKesson Brands – Aneroids and Inflation Systems - Connector Change

Overview

Due to new regulatory guidance from the FDA, all aneroids and inflation systems were changed to bayonet connectors.

Which products were affected by the Connector Change?

Item #	Description	Description 2
01-750W-11ABKGM	Wall Aneroid Sphygmomanometer	Adult, Black, McK PB
01-752M-11ABKGM	Mobile Aneroid Sphygmomanometer	Adult, Black, McK PB
01-865-10SARBGGM	Reusable Blood Pressure Cuff	Sm Adult, Royal Blue, McK PB
01-865-11ABKGM	Reusable Blood Pressure Cuff	Adult, Black, McK PB
01-865-12XBDGM	Reusable Blood Pressure Cuff	Lrg Adult, Burgundy, McK PB
01-865-13TBRGM	Reusable Blood Pressure Cuff	Thigh, Brown, McK PB
01-865-9CGRGM	Reusable Blood Pressure Cuff	Child, Green, McK PB

If customers ordered a clock aneroid and inflation system after September 2022, both arrive with the updated connector style and can be used together immediately. Nothing else is required. The connector update was implemented on 9/12/22 for these items. Customers who ordered after this time should have received the updated connector and will know by the information card provided with the product and the bayonet style connector on the product.

Note: If your coil tubing is 8974F – this will work with all new cuffs



If customers ordered new inflation systems to use with a legacy clock aneroid instrument, they’ll need to convert one of the connectors.

If the coil tubing is 891M – you will need to use 891F and can swap the cuff connector. Check the inflation system packaging for this replacement legacy connector and instructions (#891FKIT) so you can quickly retrofit the cuff in the field.

The retrofit part is included in products built up until 3/15/23. After this time, customers will need to request a retrofit kit by calling the Customer Service line for McKesson Aneroids at [1-800-232-2670](tel:1-800-232-2670).



Where can the customers get the connector if they don’t have one and need to convert?

Replacement connectors and spare parts for customers and end users who need to convert their devices, can call the customer service line for McKesson Aneroids at [1-800-232-2670](tel:1-800-232-2670).