

Medline OPA/28

High-level disinfectant solution



Work faster,
lasts longer



Safety counts and minutes matter

Time and safety are precious in processing. Getting more done in less time and reducing disinfectant exposure is easier with Medline OPA/28. It works faster and lasts longer than most OPA brands to provide safety, convenience and value.

Medline OPA/28 by the numbers

5 minutes

AER disinfection time at 25°C

2 final rinses

needed in Medivators® AER

10 minutes

Manual disinfection time at 20°C

3 final rinses

needed for manual application

28 days

reuse period

Lasts longer

Medline OPA/28 has a reuse period of 28 days—up to twice as long as other OPA brands. One bottle can do the work of two over a 28-day period, creating cost-savings and safely limiting worker exposure to disinfectants.¹



Works fast

Destroy pathogens in just 10 minutes (at room temperature) with this powerful, ready-to-use disinfectant. Some reprocessors of semi-critical medical devices use the active ingredient, 0.575% ortho-phthalaldehyde (OPA), to safely disable harmful pathogens, like:²

- MRSA
- VRE
- TB
- HIV
- Hepatitis Viruses



Bactericidal, fungicidal and tuberculocidal action in 10 minutes².




Increases cleaning cycles

Regain an extra disinfection cycle per hour with Medline OPA/28's fast cycling time². This can expand your operational efficiencies and reduce medical instrument downtime.

Other OPA Brands
(12 min cycle)

Medline OPA /28
(10 min cycle)



Get fast MRC verification

Test Strips are designed specifically for Medline OPA/28.

- **Fast**—3-second dip time; 90-second read time
- **Intuitive**—Color coding provides pass/fail reference
- **User friendly**—Large, easy-to-read indicator pad



References:
1. Provided that MRC meets requirements and product is used for full 28-day period.
2. Data available upon request.

Medline OPA/28 High-level Disinfectant

Discover the more efficient method to reprocess heat-sensitive, semi-critical medical devices that are unsuitable for sterilization.

Item Number	Description	Pkg.
MDS880PA28	High-level disinfectant, 1-gallon bottles	4/cs
MDS880PASTP	Test Strips, 50 strips/bottle	2/cs



Learn more. Contact your Medline Representative or call 1-800-MEDLINE.



**We make
healthcare
run better™**

Medline Industries, LP
Three Lakes Drive, Northfield, IL 60093
Medline United States | 1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-268-2848 | medline.ca | canada@medline.com
Medline México
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    