

McKESSON

Lumbar Support Cushion compressed

FOAM | PROVIDES LOWER BACK SUPPORT | IMPROVES SEATING POSTURE | BLUE | SINGLE PATIENT USE

For use with MFR #s 146-RTL-1493COM.


Cojín de soporte lumbar comprimido

ESPUMA | BRINDA SOPORTE PARA LA ESPALDA BAJA | MEJORA LA POSTURA AL ESTAR SENTADO | AZUL | PARA USO DE SOLO UN PACIENTE


Coussin pour soutien lombaire comprimé

MOUSSE | OFFRE UN SOUTIEN AU BAS DU DOS | AMÉLIORE LA POSTURE ASSISE | BLEU | UTILISATION UNIQUE


INSTRUCTIONS FOR USE / INSTRUCCIONES / INSTRUCTIONS



1. Unroll compressed cushion.
1. Desenrolle la almohada comprimida.
1. Dérouler le coussin comprimé.



2. Removing product from the plastic packaging will start the decompression process.
2. Retirar el producto del empaque plástico hará que comience el proceso de descompresión.
2. Le retrait du produit de l'emballage en plastique va commencer le processus de décompression.



3. Cushion will decompress to full size in 1-2 hours.
3. El cojín se descomprimirá a su tamaño máximo en 1 o 2 horas.
3. Le coussin va se décompresser jusqu'à pleine taille en 1 à 2 heures.

Not made with natural rubber latex.

Questions? Call 1-800-777-4908

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Distributed by McKesson Medical-Surgical Inc.
Richmond, VA 23233
PVN A0621
Made in China