



**WRAP YOURSELF**

**IN 360° PROTECTION.**

**Wrap-Around Reusable Isolation Gown**  
Complete front and back coverage with easy on and off design

# MORE COVERAGE, COMFORT AND SAVINGS.

You really can have it all with our Wrap-Around reusable isolation gown. Compared to disposable options, it's lighter and easier to use. Plus, it delivers exceptional savings all while providing AAMI Level 1 protection with our innovative 360-degree design.

## Feel safer and more secure.

- » **All-around protection**—complete front and back coverage
- » **Hassle-free donning and removal**—no snaps or ties
- » **Outstanding fluid resistance**—proprietary non-linting polyester
- » **Long-lasting barrier properties**—tested to last up to 75 processings
- » **More cost savings**—save up to 50% with reusable textiles over disposable gowns<sup>1</sup>
- » **Exceptional sustainability**—reduce medical waste substantially<sup>2</sup>

Item No.	Color	Size	Sweep	Length	Pkg.
MDTISO3AHY	Yellow	OSFM	72.5"	44"	dz
MDTISO3AHG	Gray				dz
MDTISO3AHYCAR	Carbon				dz
MDTISO3AHY3XL	Yellow	3XL	90.5"	44"	dz
MDTISO3AHG3XL	Gray				dz
MDTISO3AHCAR3	Carbon				dz



Don our gown three ways, see one recommended method below:



Simply slip arms into traditional sleeves



Slide left arm through opening in protective front panel

**Try one on today!**

Contact your Medline representative for more information.







Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com

FOLLOW US    

**References:** 1. Why Choose Reusable Isolation Gowns? American Reusable Textile Association. [http://www.arta1.com/sites/arta1.com/files/M7-103\\_Low-res\\_solation\\_gowns\\_sheet709.pdf](http://www.arta1.com/sites/arta1.com/files/M7-103_Low-res_solation_gowns_sheet709.pdf). Accessed August 10, 2016. 2. A comparison of reusable and disposable perioperative textiles: sustainability state-of-the-art 2012. Anesthesia and Analgesia, 2012 May;114(5):1055-66. Available at <http://www.ncbi.nlm.nih.gov/pubmed/22492184>. Accessed August 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.  
© 2016 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc.

MKT19120899 / LIT926R / 2.5M / DP / 63