



**GET DRY.**

**STAY DRY.**

**Quick Dry Reusable Underpads**  
Create the ideal skin environment


# The fast way to keep patients dry, comfortable and safe.

More than half of all older Americans struggle with incontinence.<sup>1</sup> Moisture can lead to painful and expensive skin breakdown, and even pressure ulcers.<sup>2</sup>

Introducing new **Quick Dry underpads**. The ultra-absorbent face fabric quickly pulls moisture away from the body so your patients feel dry in minutes.

- » **Rapidly absorbs moisture** to keep skin dry.
- » **Performs like a high-end disposable underpad** but with the sustainability and cost benefits of a reusable.
- » **Lays flat on the bed;** no bunching under patient.

## How Quick Dry Works



- A Face Fabric**  
Soft 100% brushed polyester with an ultra-absorbent coating to resists stains.
- B Soaker Layer**  
Heavy 10 oz. absorbency capacity. Unique fiber blend dries in minutes.
- C Barrier**  
Breathable polyurethane film helps block fluid, and allows moisture vapor to escape.
- D Base Layer**  
100% polyester blue backing anchors the pad to the bed, and lays flat.



Item No.	Description	Pkg.
MDTIU8SEFBLU	Quick Dry 100% polyester reusable underpad, 34" x 36"	2 dz/cs

**Trial QuickDry today. Contact your Medline Representative or call 1-800 MEDLINE to get started.**



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com

FOLLOW US 

1. U.S. Centers for Disease Control and Prevention's National Center for Health Statistics, WebMD. <http://www.webmd.com/urinary-incontinence-oab/news/20140625/over-half-of-seniors-plagued-by-incontinence-cdc#1>. Accessed July 12, 2017. 2. Lachenbruch, Charlie; Ribble, David; Emmons, Kirsten; VanGilder, Catherine. Pressure Ulcer Risk in the Incontinent Patient: Analysis of Incontinence and Hospital-Acquired Pressure Ulcers From the International Pressure Ulcer Prevalence™ Survey. *Journal of Wound Ostomy & Continence Nursing*; May/June 2016 - Volume 43 - Issue 3 - p 235-241. [https://journals.lww.com/jwcnonline/Abstract/2016/05000/Pressure\\_Ulcer\\_Risk\\_in\\_the\\_Incontinent\\_Patient\\_.3.aspx](https://journals.lww.com/jwcnonline/Abstract/2016/05000/Pressure_Ulcer_Risk_in_the_Incontinent_Patient_.3.aspx). Accessed July 12, 2018. Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2018 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. MKT1886223 / e18246 / PDF / 39