

# Sharps Containers

Disposal systems for every need.



**Safety**  
**at your**  
**fingertips**



# We make it easy and safe to dispose of sharps products

Each year, an estimated **600,000–800,000 needlestick injuries** and other percutaneous injuries occur among healthcare workers. The good news is the Centers for Disease Control and Prevention (CDC) estimates that 62–88 percent of these injuries can be prevented simply by using safer medical devices.<sup>1</sup>

## That's why we make our sharps containers to be among the safest in the market:

- Easy-to-see fill lines to reduce over-filling and accidental needlesticks
- Innovative lid designs offer more variety for all types of sharps and limit unwanted access to container contents
- Containers lock easily for final disposal and safe transportation

Most importantly, since we control our own supply chain, we can reduce costs and reliably deliver these products where and when you need them.

Look for this key  
throughout this  
brochure to easily  
identify the sharps  
container you need.

### Sharps container key



Red—Biohazard



Blue—Pharmaceutical



Yellow—Chemotherapy



Black—RCRA

**This table indicates general guidelines for proper waste stream segregation.**

**Regulated biohazardous medical waste**

**Hazardous waste**



**Sharps**

- Needles
- Syringes
- Pipettes
- Glass slides
- Sharp instruments
- Other single-use sharps

**Pharmaceutical waste**

- Pharmaceutical creams, ointments
- Partial vials, ampules
- Partially filled syringes
- Discontinued IVs with non-hazardous pharmaceuticals

**Chemotherapy waste**

- Trace chemotherapy waste
- Syringes needles
- Vials
- Trace chemotherapy sharps

**Hazardous RCRA**

- Hazardous compatible RCRA pharmaceutical waste
- P, U, D listed RCRA pharmaceutical waste
- Non-RCRA bulk chemo
- Hazardous aerosol pharmaceutical waste
- Pressurized Inhalers

Autoclave incineration steam sterilization and maceration (SSM) technology (red bag solutions) gringer/shredder



Low-temperature incinerator



High-temperature incinerator



Landfill



Landfill



Lined hazardous waste landfill

# Patient room sharps containers



- Designed for use in patient rooms, on med carts or in high-traffic areas
- Counter-balanced lid facilitates safe disposal
- Lid design facilitates one-handed operation
- Horizontal drop design for maximum capacity to fill line
- Easy-to-see fill levels; multiple cabinet and bracket options

Item No.	Description	Accessories	Pkg.
MDS701202F	2-gal., Translucent Red, Counter-Balance Lid	I, J	20/cs
MDS701201F	1-gal., Translucent Red, Counter-Balance Lid	I, J	32/cs
MDS705153	5-qt., Translucent Red, Counter-Balance Lid	A, C, F, G	20/cs
MDS705154	5-qt., Clear, Counter-balance Lid	A, C, F, G	20/cs
MDS705203	3-gal., Translucent Red, Counter-Balance Lid	B, C, K, L	12/cs
MDS705204	3-gal., Red Flap	D	10/cs
MDS705205	3 gal, Clear, Flap Top	D	10/cs
MDS705254	5.4-qt., Clear, Counter-Balance Lid	H	20/cs
MDS705253	5.4-qt., Red, Counter-Balance Lid	H	20/cs



MDS701201F

## Accessories

Cabinets and wall brackets free up counter space and their natural design complements any setting.

### Cabinets and glove box holder


Key	Item No.	Description	Pkg.
<b>A</b>	MDS707953	Locking Cabinet for 5-qt. Sharps Container	5/cs
<b>B</b>	MDS707930	Locking Cabinet for 3-gal. Sharps Container	2/cs
<b>C</b>	MDS707212	Glove Box Holder for 5-qt. and 3-gal. Cabinet	10/cs
<b>D</b>	MDS707940	Locking Wall Mount for Mounting MDS705204, MDS705204B, and MDS705205	5/cs



MDS705153 inside of MDS707953

### Wire wall brackets

<b>F</b>	MDS707316	Non-locking Bracket for 5-qt. Sharps Container	5/cs
<b>G</b>	MDS707315	Locking Bracket for 5-qt. Sharps Container	5/cs
<b>H</b>	MDS707317	Locking Bracket for 5.4 qt. Sharps Container	5/cs
<b>I</b>	MDS707021	Non-locking Bracket for 1- and 2-gal. Sharps Container	5/cs
<b>J</b>	MDS707112	Locking Bracket for 1- and 2-gal. Sharps Container	5/cs
<b>K</b>	MDS707003	Non-locking Bracket for 3-gal. Sharps Container	10/cs
<b>L</b>	MDS707133	Locking Bracket for 3-gal. Sharps Container	10/cs



MDS707317

# Phlebotomy sharps containers



- Fit into blood drawing trays for point-of-use disposal with a square base for stability
- Translucent top allows for easy identification of contents
- Choice of single or dual openings to accommodate a variety of sharps
- Containers lock for final disposal and safe transport
- Brackets available to mount containers to wall and free up counter space

Item No.	Description	Accessories	Pkg.
MDS705110	1-qt., Red, Phlebotomy Container	E	100/cs
MDS705115	1.5-qt., Red, Transportable Container	—	20/cs

## Accessories

Wall brackets free up counter space.

### Wire wall brackets

Key	Item No.	Description	Pkg.
<b>E</b>	MDS707011	Non-locking Bracket for 1-qt. Sharps Container	5/cs



# Multi-purpose biohazard sharps containers



- For use in restricted access areas; nestable containers save space and reduce shipping costs
- Non-torturous lid design accommodates a variety of sharps sizes
- Horizontal drop maximizes use of container volume
- Containers lock for final disposal and safe transport
- Brackets available to mount containers to wall and free up counter space

Item No.	Description	Accessories	Pkg.
MDS705201	1-gal., Red, Slide Lid	I, J	32/cs
MDS701201	1-gal., Translucent Red, Slide Lid	I, J	32/cs
MDS705202	2-gal., Red, Slide Lid	I, J	20/cs
MDS701202	2-gal., Translucent Red, Slide Lid	I, J	20/cs

## Accessories

Wall brackets free up counter space.

### Wire wall brackets

Key	Item No.	Description	Pkg.
<b>I</b>	MDS707021	Non-locking Bracket for 1- and 2-gal. Sharps Container	5/cs
<b>J</b>	MDS707112	Locking Bracket for 1- and 2-gal. Sharps Container	5/cs



# Large biohazard sharps containers




- Use in restricted access areas like operating rooms and other areas with high-waste generation
- Variety of lid designs offer new-to-market options for all sizes of sharps
- Integrated side handles make the container easy to carry
- Nestable containers save space and reduce shipping costs
- Rolling carts available for easy transport of large containers

Item No.	Description	Accessories	Pkg.
MDS705208HN	8-gal., Red, Hinge Lid	M	10/cs
MDS705208SL	8-gal., Red, Slide Lid	M, O	10/cs
MDS705208ST	8-gal., Red, Star Lid	M	10/cs
MDS705210HN	10-gal. Red, Hinge Lid	M	6/cs
MDS705210SL	10-gal., Red, Slide Lid	M, O	6/cs
MDS705210ST	10-gal., Red, Star Lid	M	6/cs
MDS705212HN	12-gal., Red, Hinge Lid	M	10/cs
MDS705212SL	12-gal., Red, Slide Lid	M, O	10/cs
MDS705212ST	12-gal., Red, Star Lid	M	10/cs
MDS705218HN	18-gal., Red, Hinge Lid	M	5/cs
MDS705218SL	18-gal., Red, Slide Lid	M, O	5/cs



Hinge lid



Slide lid



Star lid

## Accessory

### Rolling cart

Key	Item No.	Description	Pkg.
<b>M</b>	MDS707CART	Rolling Trolley Cart for all 8-18 gal. Containers	1 ea
<b>O</b>	MDS707SLIDEDECART	Rolling Foot Pedal Cart for 8-18 gal. Slide Lid Containers	1 ea

MDS705210HN

MDS705218HN

MDS705212ST

MDS705208SL



# Large waste containers



- These waste containers are designed for use in restricted access areas
- Provide safe disposal for operating room and other clinical areas with high waste generation
- Containers are nestable to help save space and reduce shipping costs
- Choice of lid designs offers options for different sizes and shapes of product to be disposed of

Item No.	Description	Accessories	Pkg.
MDS706318HN	18-gal., Yellow, Hinge Lid	M	5/cs
MDS706318SL	18-gal., Yellow, Slide Top	M, O	5/cs
MDS707318HN	18-gal., Blue, Hinge Top	M	5/cs
MDS707318SL	18-gal., Blue, Slide Top	M, O	5/cs

## Accessories

### Rolling carts

Key	Item No.	Description	Pkg.
<b>M</b>	MDS707CART	Rolling Trolley Cart for all 8–18 gal. Containers	1 ea
<b>O</b>	MDS707SLIDEDECART	Rolling Foot Pedal Cart for 8–18 gal. Slide Lid Containers	1 ea



MDS706318HN

MDS707318SL

# Large PG-II and waste containers



- PG-II containers comply with Department of Transportation (DOT) requirements to eliminate the need for secondary packaging upon disposal
- For use in restricted access areas
- Safe disposal for operating room, chemotherapy and other medical waste
- Trolley cart with foot pedal allows for hands-free operation
- Lids feature a gasket seal to resist leaking and integrated side handles for easy carrying

Item No.	Description	Accessories	Pkg.
<b>Biohazard waste</b>			
MDS705208	8-gal., Red, Flat Lid, PG-II	M, N	9/cs
MDS705308	8-gal., Red, Port Lid, PG-II	M, N	9/cs
MDS705212	12-gal., Red, Flat Lid, PG-II	M, N	8/cs
MDS705312	12-gal., Red, Port Lid, PG-II	M, N	8/cs
MDS705218	18-gal., Red, Flat Lid, PG-II	M, N	7/cs
MDS705318	18-gal., Red, Port Lid, PG-II	M, N	7/cs

<b>Chemotherapy waste</b>			
MDS705256	5.4 qt, Yellow, Counter-Balance Lid	H	20/cs
MDS706902	2-gal., Yellow, Hinge Lid	I, J	20/cs
MDS706308	8-gal., Yellow, Port Lid, PG-II	M, N	9/cs
MDS706312	12-gal., Yellow, Port Lid, PG-II	M, N	8/cs
MDS706218	18-gal., Yellow, Flat Lid, PG-II	M, N	7/cs
MDS706318	18-gal., Yellow, Port Lid, PG-II	M, N	7/cs

<b>Pharmaceutical waste</b>			
MDS705255	5.4 qt, Blue, Counter-Balance Lid	H	20/cs
MDS707902	2-gal., Blue, Hinge Lid	I, J	20/cs
MDS705204B	3-gal., Blue, Flap	D	10/cs
MDS707308	8-gal., Blue/White, Port Lid, PG-II	M, N	9/cs
MDS707318	18-gal., Blue/White, Port Lid, PG-II	M, N	7/cs

<b>RCRA hazardous waste</b>			
MDS708902	2-gal., Black, Hinge Lid	I, J	20/cs
MDS708308	8-gal., Black, Port Lid, PG-II	M, N	9/cs
MDS708318	18-gal., Black, Port Lid, PG-II	M, N	7/cs



MDS705318

MDS705208



MDS706218

MDS706308



MDS707318

MDS707308



MDS708318

MDS708308



## Accessories

### Cabinets and glove box holder

Key	Item No.	Description	Pkg.
<b>D</b>	MDS707940	Locking Wall Mount for Mounting MDS705204, MDS705204B, and MDS705205w	5/cs

### Wire wall brackets

Key	Item No.	Description	Pkg.
<b>H</b>	MDS707317	Locking Bracket for 5.4 qt. Sharps Container	5/cs
<b>I</b>	MDS707021	Non-locking Bracket for 1- and 2-gal. Sharps Container	5/cs
<b>J</b>	MDS707112	Locking Bracket for 1- and 2-gal. Sharps Container	5/cs



MDS707112

### Rolling carts

Key	Item No.	Description	Pkg.
<b>M</b>	MDS707CART	Rolling Trolley Cart for all 8–18 gal. Containers	1 ea
<b>N</b>	MDS705218CART	Rolling Foot Pedal Cart for 8–18 gal. PG-II Containers	1 ea



MDS705218HN and MDS707CART  
See page 7



MDS705212SL and MDSSSLIDECART  
See page 7



MDS705218 and MDS705218CART  
See page 9



**Reduce the risk.** For additional information, contact your Medline Representative or visit [medline.com](https://www.medline.com).

**References: 1.** Bloodborne Pathogens and Needlestick Prevention, OSHA website. Available at <https://www.osha.gov/bloodborne-pathogens/evaluating-controlling-exposure>. Accessed Dec. 7, 2022.



**We make  
healthcare  
run better™**

**Medline Industries, LP | 1-800-MEDLINE (633-5463) | [www.medline.com](https://www.medline.com)**

©2024 Medline Industries, LP. All rights reserved. Medline is a registered trademark of Medline Industries, LP. All product names, regardless of whether or not they appear with a trademark symbol, are the property of their respective companies unless otherwise noted.  
MKT19W6573392 / e22477 / 89