

McKesson Disinfectants and Cleaners/Lubricants

To learn more or place an order, please contact your McKesson Medical-Surgical Account Manager today.

Pro-Tech RTU Disinfectant Cleanser

- Effective Bactericide, Fungicide and Virucide in the presence of organic soil (5% blood serum).
- For use in pre-cleaning and early-stage decontamination of critical or semi-critical medical devices prior to sterilization or high-level disinfection.



| | QUANTITY | REORDER NO. | EC ITEM NO. | ITEM NO. |
|--------|----------|-------------|-------------|----------|
| 32 oz. | 12/cs | 53-28564 | 28644100 | 484484 |
| Gallon | 4/cs | 53-28561 | 28614100 | 484483 |

Pro-Tech Citrus Disinfectant Spray

- General purpose disinfectant/deodorant.
- Pleasant fragrance keeps air smelling fresh by counteracting objectionable odors.
- For use on non-critical medical devices prior to sterilization or high-level disinfection.
- Medical Device Uses: Tuberculocidal, Fungicidal, Bactericidal, Staphylocidal, Pseudomonacidal, Virucidal against HIV-1 (AIDS Virus), Salmonellacidal.



| | QUANTITY | REORDER NO. | EC ITEM NO. | ITEM NO. |
|--------|----------|-------------|-------------|----------|
| 16 oz. | 12/cs | 53-28584 | 28844100 | 484482 |

Multi-Enzymatic Cleanser

- Contains a synergistic combination of multiple enzymes that effectively break down blood, proteins, fat, carbohydrates and synthetic liquids, commonly found on flexible and rigid scopes, surgical instruments and other general healthcare equipment.
- Bacteriostatic, neutral pH, non-abrasive, 100% biodegradable.
- Non-foaming formula specifically created to not harm metal, plastic, rubber or corrugated tubing.
- Effective in both manual and automated reprocessing systems, as well as ultrasonic units and washer sterilizer/decontaminators.



| | QUANTITY | REORDER NO. | EC ITEM NO. | ITEM NO. |
|--------|----------|-------------|-------------|----------|
| Gallon | 4/cs | 53-28501 | 85014100 | 484478 |

Multi-Enzymatic Foaming Cleanser - RTU

- An effective pre-decontamination spray for surgical instruments and scopes.
- Thick foaming spray effectively covers the instruments and scopes, cleaning away blood, tissue, mucus and other protein-rich body fluids on contact.
- Safe for instruments and scopes.
- Neutral pH, non-abrasive, and 100% biodegradable.



| | QUANTITY | REORDER NO. | ITEM NO. |
|--------|-----------------|--------------------|-----------------|
| 32 oz. | 12/cs | 53-28517 | 484479 |

Low Suds Liquid Instrument Detergent

- Designed to effectively clean medical devices and surgical instruments.
- Low-foaming, neutral pH formula is appropriate for use in automated washers, as well as manual cleaning systems.
- Non-corrosive, contains no chlorines, phosphates or hydroxides and will not harm surgical instruments.



| | QUANTITY | REORDER NO. | ITEM NO. |
|--------|-----------------|--------------------|-----------------|
| Gallon | 4/cs | 53-28551 | 484481 |

Instrument Lubricant and Rust Inhibitor

- Effective protection for instrument surfaces, joints, and cutting edges.
- Non-silicone formula.
- Does not interfere with sterilization process.
- Coats metal instruments with an oil-based film that lubricates the instrument while also inhibiting corrosion.
- Helps prevent instrument spotting, rusting and staining and helps reduce instrument freeze-up by keeping hinges working freely.



| | QUANTITY | REORDER NO. | ITEM NO. |
|--------|-----------------|--------------------|-----------------|
| Gallon | 4/cs | 53-28571 | 484480 |

The product information contained in this flier, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers, not McKesson Medical-Surgical or its affiliates ("McKesson") and have not been independently verified by McKesson. McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the flier, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print fliers by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this flier is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical

8741 Landmark Road
Richmond, VA 23228

www.mckesson.com

800.446.3008

©2008 McKesson Medical-Surgical Inc.
2008-0086