

Product Change Notification**Date of Publication: October 10, 2016****Characteristic of Change: New Purple Chemistry**

Valued Customer,

We are pleased to announce an improved formulation of McKesson OPA/28 High-Level Disinfectant. You will now notice that the solution is tinted purple, for easier identification and differentiation from other clear liquids, including water. You may also notice a change in the foaming characteristics of McKesson OPA/28 depending on your specific use. We have updated our packaging to indicate that the solution is now tinted purple.

These improvements are well-tested to optimize the user experience while ensuring the same efficacy, speed, and re-use life that you've come to expect. Instructions for use, label claims, and key performance indicators including concentration of OPA and operating pH are identical to the original formulation. The clear-colored McKesson OPA/28 from previous lots can be mixed with or used to top off the purple-tinted McKesson OPA/28 without effect on the efficacy of the product.

We thank you for being a loyal customer and are excited to be able to further enhance McKesson OPA/28. If you have any questions about McKesson OPA/28 High-Level Disinfectant, please contact your local McKesson Medical-Surgical Inc. sales representative.

Distributed by McKesson Medical-Surgical Inc.
9954 Mayland Drive
Richmond, VA 23233

www.mckesson.com