




## **Flex-Grip™** Cane Tip

### Instructions:

1. Remove the old tip from your cane.
2. Slip on the Flex-Grip™ cane tip, making sure that the Airgo logo is aligned with the front of the cane (the side that is facing forward when you walk). This will ensure optimal grip and traction for a more natural walking motion.

### Fits most aluminum canes

To verify that the Airgo® Flex-Grip™ cane tip is right for your cane, remove the old cane tip and match the shaft end with the circle below.



Questions? Comments?  
[www.airgomobility.com](http://www.airgomobility.com)

Distributed by / Distribué par  
Drive Medical Canada  
256 Aviva Park Drive  
Vaughan, Ontario Canada L4L 9C7  
[www.drivemedical.com](http://www.drivemedical.com)  
USA: 866-694-5085  
CAN: 866-752-1363

REV4.12.07.16 P/N 921-735001

### SPECIFICATIONS

Project: **Airgo Flex-Grip Instructions Rev3 0815 FRONT**

Item No.: **735-001** Part No.: **921-735001**

Flat Size: **4.5" x 6.5"**

Assembled size: **4" X 6"**

Colours:  B/W  CMYK  Spots:

Special Instructions: **Double-sided Instruction Sheet**

Generated by:



Pour vous et vos besoins...



## Embout **Flex-Grip**<sup>MC</sup>

### Instructions:

1. Retirer le vieil embout de votre canne.
2. Placer l'embout Flex-Grip<sup>MC</sup> en prenant soin d'aligner le logo Airgo à l'avant de la canne, de sorte qu'il pointe devant vous lorsque vous marchez. Cela vous permettra de profiter de l'adhérence accrue que vous offre l'embout Flex-Grip<sup>MC</sup> et, ainsi, de marcher plus confortablement.

### S'ajuste à la plupart des cannes en aluminium

Pour vérifier que l'embout Flex-Grip<sup>MC</sup> convient à votre canne, retirez le vieil embout puis placez le bout de votre canne sur le cercle gris ci-dessous.



Questions? Commentaires?  
[www.airgomobility.com](http://www.airgomobility.com)

Distributed by / Distribué par  
Drive Medical Canada  
256 Aviva Park Drive  
Vaughan, Ontario Canada L4L 9C7  
[www.drivemedical.com](http://www.drivemedical.com)  
USA: 866-694-5085  
CAN: 866-752-1363

### SPECIFICATIONS

Project: **Airgo Flex-Grip Instructions Rev3 0815 BACK**

Item No.: **735-001** Part No.: **921-735001**

Flat Size: **4.5" x 6.5"**

Assembled size: **4" X 6"**

Colours:  B/W  CMYK  Spots:

Special Instructions: **Double-sided Instruction Sheet**

Generated by: