

Post-mortem bags and shrouds

We meet your needs for assortment, value and quality.



Body bags*

- Intended for containment use only, not designed for transportation
- Heat sealed edges allow for leak resistance and strength
- Not weight-bearing, must be supported during lifting/repositioning
- Ideal for morgues, hospitals and EMS services

Item No.	Description	Thickness	Maximum Weight	Pkg.
NON70540WM	Adult, 36"W x 90"H, PVC, c-curve zipper, White	8 mil	200 lbs.	10/cs
NON70548WM	X-large, 48"W x 90"H, PVC, c-curve zipper, White	8 mil	200 lbs.	10/cs
NON70572WM	X-large, 72"W x 96"H, PVC, c-curve zipper, White	10 mil	450 lbs.	5/cs
NON70550WM	Child, 28"W x 48"H, PVC, straight zipper, White	8 mil	150 lbs.	30/cs
NON70560WM	Infant, 18"W x 28"H, PVC, straight zipper, White	8 mil	100 lbs.	30/cs



NON70540WM

Transport bags*

- Bags are reinforced with six fabric handles to assist in repositioning and lifting
- Heat-sealed finish prevents leakage
- Use disaster bag (see page 2) if larger weight capacity is needed

Item No.	Description	Thickness	Maximum Weight	Pkg.
NON70740W	Adult, 36"W x 90"H, PVC, c-curve zipper, White	15 mil	500 lbs.	5/cs



NON70740W

*All bags and kits include toe tags

Post-mortem bags continued

Disaster bags*

- Thickest and most durable bags we offer
- Ten sewn-in, ergonomic handles provide ease in handling and stability in the most demanding environments
- Heavy-duty thickness maximizes weight-bearing capabilities

Item No.	Description	Thickness	Maximum Weight	Pkg.
NON80540W	Adult, 36"W x 96"H, PVC, c-curve zipper, Black	22 mil	650 lbs.	5/cs
NON80549W	X-large, 76"W x 96"H, PVC, c-curve zipper, Black	22 mil	850 lbs.	2/cs



NON80540W

PVC/DEHP-free bags*

- For facilities who want PVC/DEHP-free alternatives
- PE scrim bag fulfills chlorine-free burn requirements
- PE scrim and DEHP-free bags are phthalate free

Item No.	Description	Thickness	Maximum Weight	Pkg.
NON60540B	40"W x 91"H, PE scrim, c-curve zipper, Black	5.5 mil	200 lbs.	10/cs
NON60540W	36"W x 90"H, DEHP-free, c-curve zipper, White	8 mil	200 lbs.	10/cs



NON60540W

*All bags and kits include toe tags

Shroud packs and kits

Shroud packs*

- Accessories to prepare the cadaver for final arrangements and provide coverage from view
- **Adult pack includes:** one 54"W x 108"L plastic sheet, three 36" ties, two 60" ties, three ID tags, two absorbent pads and one chin strap
- **X-large pack includes:** one 72"W x 108"L plastic sheet, three 36" ties, two 60" ties, three ID tags, two absorbent pads and one chin strap

Item No.	Description	Pkg.
NON70500	Adult, 54"W x 108"L	24/cs
NON70510	X-large, 72"W x 108"L	24/cs



Shroud kits*

- Maximize efficiency and convenience with our kits
- Essential for standardizing product usage and adherence to facility protocols
- Kits include all components listed in our shroud packs, plus a body bag
- Not weight-bearing; must be supported during lifting/repositioning

Item No.	Description	Thickness	Maximum Weight	Pkg.
NON90540WM	Adult, 36" W x 90" H, c-curve zipper	8 mil	200 lbs.	10/cs
NON90549WM	X-large, 72" W x 96" H, c-curve zipper	8 mil	200 lbs.	5/cs
NON70911WM	Adult, 36" W x 90" H, wrap-around zipper	10 mil	450 lbs.	10/cs




NON90540WM

*All bags and kits include toe tags

Shroud product information and guidelines


- Reduces the risk of bodily fluid leakage
- Allows the body to remain sturdy during repositioning and transportation
- Provides total coverage from view




1. Set the cadaver on the plastic shroud sheet and place the absorbent pads under the rear end




2. Fasten chin strap and tie at the top of the head




3. Fold arms over abdomen and couple wrists and ankles together; gently use the ties to secure



4. Attach the identification tag to the toe



5. Wrap the cadaver in the plastic shroud sheet and use the remaining three ties to tie above the elbows, at waist, and below knees; fasten ID tag on a tie



6. Place cadaver inside the body bag, secure shut, and tie the remaining ID tag to the outer zipper

For samples, pricing and availability
contact your Medline Representative.