

Alcohol Prep Pads

McKesson Alcohol Prep Pads help to prepare the skin prior to injection or venipuncture. They are composed of a soft, 2-ply material designed for the administration of blood tests and insulin or heparin injections. McKesson Alcohol Prep Pads are available in both regular and lightweight, as well as in a variety of sizes, for caregiver convenience. Regular alcohol prep pads are constructed of a non-woven material, while the lightweight are of a lesser weight to minimize waste.

Features and Benefits

- First aid antiseptic designed for preparation of the skin prior to injection or venipuncture
- Enclosed in foil pouches to maintain sterility
- Impregnated with 70% isopropyl alcohol
- Not made with natural rubber latex



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Alcohol Prep Pads, Medium, 1.2" x 2.6"	58-204	200ea/bx, 20bx/cs
Alcohol Prep Pads, Large, 1.8" x 3.5"	58-404	100ea/bx, 10bx/cs
Lightweight Alcohol Prep Pads, Medium, 1.2" x 2.6"	16-527	200ea/bx, 20bx/cs
Lightweight Alcohol Prep Pads, Large, 1.8" x 3.3"	16-528	100ea/bx, 10bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.