

## Urgent Assist Bags

McKesson urgent assist bags allow for easy, hygienic disposal of vomit or urine while eliminating messes and clean-up. Caregivers can easily twist the bags closed and hook them into the notched rings to help reduce odors. The twist and hook enclosure also provides a way to help prevent contamination. Our emesis bags are graduated in 4 fl oz. increments to allow for monitoring of content volume. Medical staff can hang the bag near patients to avoid any spills.

### Features and benefits

- For easy, hygienic collection and disposal of urine or vomit
- 7-inch wide rigid collar allows for easy use
- To use, press ends of collar bottom to open funnel and position over mouth or genital area
- Inner valve automatically seals contents in, preventing leaks and spills
- For added protection against spills, the funnel can be shut by folding in on its scored line
- No clean-up, easily and discreetly disposed of when convenient
- Opaque (white) bag, which conceals inner contents
- Not made with natural rubber latex



### ORDERING INFORMATION

MFR #	DESCRIPTION	ATTRIBUTES	SIZE	QTY
16-7000	Urgent Assist Bag	Opaque	40 fl oz	1 ea, 12 ea/pack, 20 packs/cs
16-7000HP	Urgent Assist Bag, Protective Hand Sleeve	Opaque	40 fl oz	1 ea, 12 ea/pack, 20 packs/cs
16-8000	Urgent Assist Bag	Graduated, Clear	40 fl oz	1 ea, 12 ea/pack, 20 packs/cs
16-8000HP	Urgent Assist Bag, Protective Hand Sleeve	Graduated, Clear	40 fl oz	1 ea, 12 ea/pack, 20 packs/cs

#### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.