

Disposable Germicidal Surface Wipes

McKESSON

McKesson disposable germicidal surface wipes are an ideal cleaning and disinfecting solution for hard, non-porous surfaces. They are made of non-woven cloth and contain a stable, low-pH-formulated disinfectant and deodorant. The wipes are pre-moistened and disposable. A self-threading dispenser is built into the lid to help ensure hassle-free dispensing. Active ingredients: N-alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chloride and N-alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chloride, and isopropanol (21%).

Features and Benefits

KILLS DANGEROUS PATHOGENS

EPA-registered to kill many viruses, bacteria, and fungi in 2 minutes on pre-cleaned, hard, non-porous surfaces at room temperature

LOW-pH DISINFECTANT FORMULA

Contains a stable, low-pH-formulated disinfectant and deodorant for effective cleaning performance

DURABLE NON-WOVEN MATERIAL

Constructed from durable, non-woven cloth for reliable strength during use

NON-ABRASIVE SURFACE COMPATIBILITY

Safe for use on nonporous surfaces without scratching or damaging finishes


CONVENIENT DISPENSING DESIGN

Packaged for easy, efficient dispensing to support workflow needs

ECONOMICAL SINGLE-WIPE DISPENSING

Self-threading dispenser lid helps ensure hassle-free, controlled, single-wipe dispensing while reducing waste

NOT MADE WITH NATURAL RUBBER LATEX



DO NOT FLUSH

Ordering Information

E1 #	MFR #	DESCRIPTION	QTY
1244103	50-66165	Disposable Germicidal Surface Wipes, Flat Pack, 7 in x 9 in	80 wipes/pk, 12 pk/cs
880563	50-66160	Disposable Germicidal Surface Wipes, Canister, 6 in x 6.75 in	160 wipes, 12/cs
949454	50-66170	Disposable Germicidal Surface Wipes, Canister, 9 in x 12 in	65 wipes, 12/cs
624463	M59963	Wall Bracket	1 ea

Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson and the McKesson logo are registered trademarks of McKesson Corporation.

© 2026 McKesson Medical-Surgical Inc.

Rev. 00 03/26

Wipes: hospital and general disinfection or bactericidal

This product kills the following bacteria in two minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

- *Acinetobacter baumannii* {ATCC 19606}
- *Escherichia coli* – with [(extended beta-lactamase resistance) [(ESBL)]] [(E. coli)] {ATCC BAA-196}
- *Enterobacter aerogenes* {ATCC 13048}
- Carbapenem-resistant *Klebsiella pneumoniae* [(CRKP)] {ATCC BAA-1705}
- *Listeria monocytogenes* {ATCC 19117}
- Methicillin-resistant *Staphylococcus aureus* [(MRSA)] {ATCC 33591}
- *Mycobacterium bovis*, BCG [(tuberculosis)][(TB)] {Organon Teknika Corp}
- *Pseudomonas aeruginosa* {ATCC 15442}
- *Salmonella enterica* {ATCC 10708}
- *Serratia marcescens* {ATCC 14756}
- *Staphylococcus aureus* {ATCC 6538}
- Vancomycin – Resistant *Enterococcus faecium* [(VRE)]{ATCC 51559}

Wipes: virucidal*

This product kills the following viruses* in two minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

- SARS-related coronavirus 2 SARS-CoV-2
- Avian influenza A reassortant virus [(H3N2)] {ATCC VR-2072}
- Rotavirus [Strain WA]
- Hepatitis B virus* [(HBV)] [Strain LeGarth]
- Hepatitis C virus* [(HCV)]
- Herpes simplex II virus [(HSV 2)] [G Strain]
- Human immunodeficiency virus* [(HIV-1)] [(associated with AIDS)]
- Influenza A virus [(H1N1)] [Strain A/Swine/Iowa/15/30] [(formerly called “swine flu”)] {ATCC VR-333}
- Influenza A virus [Strain A2/Hong Kong] {ATCC VR-544}
- Rhinovirus [Type 14] [RV-14] {ATCC VR-284}

Wipes: fungicidal*

This product kills the following fungi in two minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

- *Trichophyton interdigitale* [(formerly known as *Trichophyton mentagrophytes*)] [(ringworm)] [(Athletes Foot)] {ATCC 9533}