

## Chemotherapy Procedure Gowns

McKesson chemotherapy gowns, available in open back and full back designs, provide healthcare professionals with reliable protection while preparing or handling chemotherapy drugs. Gown materials and seams are tested against 12 chemotherapy drugs per ASTM F739-12. Our ultrasonically-sealed full back gowns are AAMI Level 2 rated for barrier performance. Choose McKesson for a wide variety of high-quality gowns and related supplies.

### Features and Benefits

- Polycoated material delivers durability and protection; soft, spunbound polypropylene is coated with a layer of polyethylene film
- 2 back designs: open back or full back
- Neck closure options: hook-and-loop, over-the-head & tie
- Secured with ties at waist for added comfort
- Low-linting material helps to provide a clean environment
- Long sleeves with knit cuffs provide a comfortable fit
- Full back gowns are rated AAMI 2
- Not made with natural rubber latex



#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (MINUTES)	ADDITIONAL OBSERVATIONS
Carmustine 3.3 mg/ml (3,300 ppm)	78.3	Slight swelling and no degradation
Cisplatin 1.0 mg/ml (1,000 ppm)	>240	Slight swelling and no degradation
Cyclophosphamide 20.0 mg/ml (20,000 ppm)	>240	Slight swelling and no degradation
Dacarbazine 10.0 mg/ml (10,000 ppm)	>240	Slight swelling and no degradation
Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm)	>240	Slight swelling and no degradation
Etoposide 20.0 mg/ml (20,000 ppm)	>240	Slight swelling and no degradation
Fluorouracil 50.0 mg/ml (50,000 ppm)	>240	Slight swelling and no degradation
Methotrexate 25.0 mg/ml (25,000 ppm)	>240	Slight swelling and no degradation
Mitomycin C 0.5 mg/ml (500 ppm)	>240	Slight swelling and no degradation
Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm)	>240	Slight swelling and no degradation
Thiotepa 10 mg/ml (10,000 ppm)	>240	Slight swelling and no degradation
Vincristine Sulfate 1.0 mg/ml (1,000 ppm)	>240	Slight swelling and no degradation

## ORDERING INFORMATION

CHEMOTHERAPY PROCEDURE GOWNS					
MFR #	E1 #	DESCRIPTION	ATTRIBUTES	SIZE	QTY
16-58KVM	1177922	Open-Back, Hook-And-Loop Neck Closure	Fluid Resistant, Set-in Sleeves, Knit Cuffs, Blue	Medium, 41 in	10/bg, 3 bg/cs
16-58KVL	1177923	Open-Back, Hook-And-Loop Neck Closure	Fluid Resistant, Set-in Sleeves, Knit Cuffs, Blue	Large, 44 in	10/bg, 3 bg/cs
16-58KVX	1177924	Open-Back, Hook-And-Loop Neck Closure	Fluid Resistant, Set-in Sleeves, Knit Cuffs, Blue	X-Large, 51 in	10/bg, 3 bg/cs
16-58KV2X	1177925	Open-Back, Hook-And-Loop Neck Closure	Fluid Resistant, Set-in Sleeves, Knit Cuffs, Blue	XX-Large, 53 in	10/bg, 3 bg/cs
16-58KOL	1177926	Open-Back, Over-The-Head	Fluid Resistant, Set-in Sleeves, Knit Cuffs, Blue	Large, 44 in	10/bg, 3 bg/cs
16-54KVM	1177927	Full-Back, Hook-And-Loop Neck Closure, AAMI Level 2	Set-in Sleeves, Knit Cuffs, Blue	Medium, 41 in	10/bg, 3 bg/cs
16-54KVL	1177928	Full-Back, Hook-And-Loop Neck Closure, AAMI Level 2	Set-in Sleeves, Knit Cuffs, Blue	Large, 44 in	10/bg, 3 bg/cs
16-55KVX	1177929	Full-Back, Hook-And-Loop Neck Closure, AAMI Level 2	Set-in Sleeves, Knit Cuffs, Blue	X-Large, 51 in	10/bg, 3 bg/cs
16-55KV2X	1177930	Full-Back, Hook-And-Loop Neck Closure, AAMI Level 2	Set-in Sleeves, Knit Cuffs, Blue	XX-Large, 53 in	10/bg, 3 bg/cs

McKesson Medical-Surgical  
 9954 Mayland Drive  
 Richmond, VA 23233  
[mms.mckesson.com/mckesson-brands](http://mms.mckesson.com/mckesson-brands)

Post-Acute Care and Home Health  
 800.347.2456

Home Medical Equipment  
 888.822.8111

Primary Care  
 866.625.2679

McKesson Brands Support  
 800.777.4908