

Drugs of Abuse

McKesson drugs of abuse test cup is a CLIA-waived urine drug screen for multiple drugs at once. This sensitive and specific test provides quick results in only 5 minutes. They come with built-in specimen validity test and temperature strip to help detect specimen tampering. Each box includes a drugs of abuse test cup procedure card, drug screen cup instruction insert, important technical bulletin, specimen validity tests procedures for cup insert, security seal stickers and adulteration color comparison chart.

Features and Benefits

MFR# 16-5125A3

- 12-drug panel with adulterants
- Test includes AMP, BAR, BZO, COC, mAMP/MET, MDMA, MOP300, MTD, OXY, PCP, TCA, THC (OX, pH, SG)
- CLIA waived

MFR# 16-6125A3

- 12-drug panel with adulterants
- Test includes AMP, BAR, BUP, BZO, COC, mAMP/MET, MDMA, MOP300, MTD, OXY, PCP, THC (OX, pH, SG)
- CLIA waived



MFR# 16-5125A3

Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Features and Benefits

MFR# 16-1145A3

- 14-drug panel with adulterants
- Test includes AMP, BAR, BUP, BZO, COC, mAMP/MET, MDMA, MOP300, MTD, OXY, PCP, PPX, TCA, THC (OX, pH, SG)
- CLIA waived

MFR# 16-8105A3

- 10-drug panel with adulterants
- Test includes AMP, BUP, BZO, COC, mAMP/MET, MDMA, MTD, MOP300, OXY, THC (OX, pH, SG)
- CLIA waived



MFR# 16-1145A3

ORDERING INFORMATION

| DRUGS OF ABUSE TEST CUPS | | |
|--------------------------|--------------------------------|-------------------|
| MFR # | DRUG PANEL | QTY |
| 16-5125A3 | 12-drug panel with adulterants | 1/ea, 25/bx, 4/cs |
| 16-6125A3 | 12-drug panel with adulterants | 1/ea, 25/bx, 4/cs |
| 16-1145A3 | 14-drug panel with adulterants | 1/ea, 25/bx, 4/cs |
| 16-8105A3 | 10-drug panel with adulterants | 1/ea, 25/bx, 4/cs |

Urine Screening Device Cutoff Levels

| TEST | CALIBRATOR | CUTOFF (ng/mL) |
|--|--|-------------------------|
| Amphetamine (AMP) Amphetamine (AMP 500) | d-Amphetamine d-Amphetamine | 1000 ng/mL 500 ng/mL |
| Barbiturates (BAR) | Secobarbital | 300 ng/mL |
| Benzodiazepines (BZO) | Oxazepam | 300 ng/mL |
| Buprenorphine (BUP, BUPG) | Buprenorphine | 10 ng/mL |
| Cocaine (COC) Cocaine (COC 150) | Benzoyllecgonine Benzoyllecgonine | 300 ng/mL 150 ng/mL |
| Marijuana (THC) | 11-nor- Δ^9 -THC-9-COOH | 50 ng/mL |
| Methadone (MTD) | Methadone | 300 ng/mL |
| Methamphetamine (mAMP, MET) Methamphetamine (mAMP 500, MET 500) | d-Methamphetamine d-Methamphetamine | 1000 ng/mL 500 ng/mL |
| Methylenedioxymethamphetamine (MDMA), Ecstasy | d,l Methylenedioxymethamphetamine | 500 ng/mL |
| Opiate (MOP, OPI 300) Opiate (OPI) | Morphine Morphine | 300 ng/mL 2000 ng/mL |
| Oxycodone (OXY) | Oxycodone | 100 ng/mL |
| Phencyclidine (PCP) | Phencyclidine | 25 ng/mL |
| Propoxyphene (PPX) | Propoxyphene | 300 ng/mL |
| Tricyclic antidepressants (TCA) | Nortriptyline | 1000 ng/mL |

Urine drug screen assays provide only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) or liquid chromatography-tandem mass spectrometry (LC-MS/MS) are the preferred confirmatory methods. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are indicated.

IMPORTANT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908