

Raised Toilet Seat

According to the Centers for Disease Control and Prevention, over 14 percent of bathroom-related injuries in older adults occur “when standing up from, sitting down on, or using the toilet.” The added height of the seat and the durable arm rails allow users to safely sit and stand back up from the toilet. McKesson's raised toilet seat enables users to avoid falls and achieve greater independence and confidence. The locking device is engineered to fit securely on almost any toilet.

Features and Benefits

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Easily attaches to toilet bowl, tool-free with front-clamp, worm-screw locking mechanism engineered to prevent toilet seat from shifting
- Heavy-duty, molded-plastic construction allows the product to remain lightweight without sacrificing strength and durability
- Foam-padded, metal arms sit 5-3/4 inches above the seat to provide security, comfort, and ease of mobility
- Packaged for retail sale
- May be eligible for Flexible Spending Account (FSA) Reimbursement
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	HEIGHT	WEIGHT CAPACITY	HCPCS	QTY
146-RTL12027RA	Raised Toilet Seats, Removable Padded Arms, Front-Locking	5 in	300 lbs	E0244	1 ea
146-RTL12064	Raised Toilet Seats	4 in	400 lbs	E0244	1 ea

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates (“McKesson”). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908