

## RSV ANTIGEN TEST

The Consult RSV antigen test is intended for the rapid, qualitative detection of respiratory syncytial virus fusion protein directly from nasopharyngeal swab, and nasal aspirate specimens in children less than 6 years and adults over the age of 60. The test is intended for use as an aid in the rapid laboratory diagnosis of acute respiratory syncytial virus infection in patients with symptoms consistent with RSV infection. It is recommended that negative test results be confirmed by cell culture.

### FEATURES AND BENEFITS

- CLIA waived
- Rapid detection of RSV in patients <6 years old or >60 years old
- Built in control features in cassette (Liquid positive and negative controls provided)
- Results in 15 minutes
- Can use nasopharyngeal swabs or nasal aspirates for testing
- For nasopharyngeal swabs - (Sensitivity = 95.8%, Specificity= 98.7%)
- For nasal aspirates - (Sensitivity = 94.6%, Specificity= 96.2%)
- Not made with natural rubber latex



E1#	MFR #	DESCRIPTION	CLIA	CONTENTS	QTY
<a href="#">1270178</a>	RSV002PBS	RSV Antigen Test	CLIA Waived	Test Cassettes, Nasopharyngeal Swabs, Extraction Tubes, Disposable Pipettes, Quick Reference Instruction, Instructions for Use, Positive and Negative Controls, Extraction Reagent, PBS Tubes	20 ea/kt 20 kt/cs

**Satisfaction Guaranteed**

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.