

## Barrier island dressings

McKesson barrier island dressings combine a non-adherent absorbent gauze pad, a water-resistant barrier layer, and a non-woven adhesive border all in one. They can be used as either a primary or secondary dressing. McKesson barrier island dressings provide an effective water barrier to permit showering/bathing.

### Features and benefits

- Offers the convenience of an all-in-one dressing by combining an absorbent non-adherent central pad surrounded by a breathable adhesive border, with an addition of a polyurethane film barrier
- The non-adherent pad in the center is designed to avoid disrupting the wound bed upon removal
- Effective water barrier to permit showering/bathing
- Primary or secondary bordered non-adherent gauze composite dressing with film backing that protects from outside contamination
- Wear time of 3-5 days or as directed by a healthcare professional
- Not made with natural rubber latex



### ORDERING INFORMATION

IMAGE	E1	MFR #	DESCRIPTION	DRESSING SIZE	PAD SIZE	QTY
	<a href="#">488921</a>	16-89166	Barrier Island Dressing Water-Resistant, Sterile	6 in x 6 in (15.2 cm x 15.2 cm)	4 in x 4 in (10.2 cm x 10.2 cm)	1/ea 25/bx 4bx/cs
	<a href="#">488920</a>	16-89144	Barrier Island Dressing Water-Resistant, Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	2 in x 2 in (5.1 cm x 5.1 cm)	1/ea 25/bx 4bx/cs

### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.