

## IV Catheters

McKesson IV Catheters are designed to be inserted into the patient's vascular system for short term use (less than 30 days) to sample blood, monitor blood pressure, or administer fluids intravenously. The thin wall needle helps assure easy puncture, while the design of the needle hub allows the caregiver to maintain catheter stability after insertion, reducing risk of blood vessel and tissue damage.

### Features and Benefits

- IV catheters are made of soft, flexible material
- Sterile, single-use and non-winged
- Designed to assist clinicians in delivering medical care to their patients
- The translucent hub enables rapid detection of blood flashback
- Color coded for easy size identification:
  - 20G Pink
  - 22G Blue
  - 24G Yellow
- Not made with natural rubber latex



### ORDERING INFORMATION

MFR #	DESCRIPTION	SPECS	QTY
6048	IV Catheter, Sterile, 22G x 1 in	22G x 1 in	1 ea, 50 ea/bx, 20 bx/cs
6049	IV Catheter, Sterile, 24G x 0.75 in	24G x 0.75 in	1 ea, 50 ea/bx, 20 bx/cs
6050	IV Catheter, Sterile, 20G x 1 in	20G x 1 in	1 ea, 50 ea/bx, 20 bx/cs
6051	IV Catheter, Sterile, 20G x 1.25 in	20G x 1.25 in	1 ea, 50 ea/bx, 20 bx/cs

#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical  
9954 Mayland Drive  
Richmond, VA 23233  
[mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands)

Post-Acute Care and Home Health  
800.347.2456

Home Medical Equipment  
888.822.8111

Primary Care  
866.625.2679

McKesson Brands Support  
800.777.4908