

Extension feeding sets

McKesson extension feeding sets connect a button-style gastrostomy tube to a pump set to allow the administration of nutrition, hydration, and medication. They are used as an extension set to other gastric feeding devices, incorporating safety connectors which may reduce the risk of accidental connection of an IV system to the enteral system, or the enteral system to the IV system. Feeding sets are available in two lengths; 12-inch or 24-inch, with two different types of secure lock connectors and one or two ENFit® feeding ports.

Features and benefits

- Single ENFit port with straight or right angle connector options for bolus feeding, flushing and medication deliver
- Dual ENFit Ports with right angle connector allows for continuous feeding, easy flushing and administration of liquid medication without disconnecting set
- Length options: 12-inch or 24-inch
- G-Tube connector options: right angle or straight
- ISO 80369-3 compliant ENFit connector designed to reduce the risk of tubing misconnections
- Compatible only with MIC Key® and MINI One® low profile gastrostomy feeding tubes
- Non-sterile
- Not made with natural rubber latex
- Not made with DEHP (diethylhexyl-phthalate), BPA (bisphenol A) or phthalate plasticizers



Extension set connection guide video:



Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

MFR #	DESCRIPTION	LENGTH	CONNECTOR	FEEDING USAGE	QTY
194-0124-12	Extension Feeding Set, Single ENFit® Port with Clamp, Non-sterile, Single Patient Use	12-inch	Right Angle	Bolus	1 bag, 5 bags/bx, 6 bx/cs
194-0124-24	Extension Feeding Set, Single ENFit® Port with Clamp, Non-sterile, Single Patient Use	24-inch	Right Angle	Bolus	1 bag, 5 bags/bx, 6 bx/cs
194-0123-12	Extension Feeding Set, Single ENFit® Port with Clamp, Non-sterile, Single Patient Use	12-inch	Straight	Bolus	1 bag, 5 bags/bx, 6 bx/cs
194-0123-24	Extension Feeding Set, Single ENFit® Port with Clamp, Non-sterile, Single Patient Use	24-inch	Straight	Bolus	1 bag, 5 bags/bx, 6 bx/cs
194-0121-12	Extension Feeding Set, Dual ENFit® Port with Clamp, Non-sterile, Single Patient Use	12-inch	Right Angle	Continuous	1 bag, 5 bags/bx, 6 bx/cs
194-0121-24	Extension Feeding Set, Dual ENFit® Port with Clamp, Non-sterile, Single Patient Use	24-inch	Right Angle	Continuous	1 bag, 5 bags/bx, 6 bx/cs

