

Fluid Solidifier

McKesson fluid solidifier provides fast and effective infectious waste fluid containment by quickly solidifying biohazardous fluids in suction canisters, emesis basins, liners and bags, kick buckets, bed pans and spills. Comprised of fast-acting superabsorbent polymers, every ounce of this product solidifies 1cc of fluid. McKesson fluid solidifier helps eliminate spilling, splashing, and aerosolization and minimizes the risk of cross contamination. This product is ideal for use in the operating environment where the majority of suction canister waste is collected.

Features and Benefits

- Provides fast and effective infectious waste fluid containment by quickly solidifying biohazardous fluids
- Screw-cap bottle for easy-pouring and secure storage
- Eliminates the need to transport biohazardous fluids in liquid form
- Reduces exposure to infectious waste
- Available in multiple configurations
- Chlorine-free formula
- Protects healthcare workers and promotes workplace safety
- Not to be used for disinfection or sanitization
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

MFR #	DESCRIPTION	GEL INFO	NET WEIGHT	QTY
547-750	Fluid Solidifier	Gels up to 750 cc	0.74 oz	1 ea, 75 ea/cs
547-1500	Fluid Solidifier	Gels up to 1,500 cc	1.48 oz	1 ea, 75 ea/cs
547-3000	Fluid Solidifier	Gels up to 3,000 cc	3.1 oz	1 ea, 75 ea/cs
547-16000	Fluid Solidifier	Gels up to 16,000 cc	16.93 oz	1 ea, 20 ea/cs
547-750P	Fluid Solidifier, Drop-In Pack	Gels up to 750 cc	0.74	1 ea, 100 ea/cs
547-1500P	Fluid Solidifier, Drop-In Pack	Gels up to 1,500 cc	1.48 oz	1 ea, 100 ea/cs

McKesson Medical-Surgical
 9954 Mayland Drive
 Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
 800.347.2456

Home Medical Equipment
 888.822.8111

Primary Care
 866.625.2679

McKesson Brands Support
 800.777.4908