

Capillary Products

McKesson Capillary Products are designed to help with the collection and organization of blood samples in healthcare facilities. The plastic micro-hematocrit capillary tubes help to eliminate the risk of injury and contamination. They have color-coded ends to help medical staff quickly and easily identify plain and heparin treated tubes. Our line includes needle holders and capillary tube holding trays to help provide a convenient method of handling needles and capillary tubes during centrifugation and disposal.

Features and Benefits

- Holding trays with 24 numbered places for easy identification and wax sealant to withstand centrifugation
- Plastic capillary tubes are non-sterile and available untreated or heparin treated
- Ammonium heparinized capillary tubes to help minimize blood clotting
- Disposable, multi-sample needle holders available as universal fit
- Not made with natural rubber latex



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Capillary Tube Holding Trays	553	10ea/bx
Micro-Hematocrit Capillary Tubes, Plastic, Untreated	554	100ea/vial, 10vials/bx
Micro-Hematocrit Capillary Tubes, Plastic, Ammonium Heparinized	555	100ea/vial, 10vials/bx
Needle Holders	556	100ea/bg, 10bg/bx

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.