

**MCKESSON**

**PREVENT<sup>®</sup>**

## Safety Pen Needles

McKesson's Prevent<sup>®</sup> safety pen needles help avoid accidental needle stick injuries from both ends of the device with McKesson's uniquely designed hub and needle shield. The user end is contained within a needle shield and the non-patient end has an elongated hub to help prevent needle contact. Works with most leading pen injectors. Our safety pen needles reduce total injection costs while minimizing patient discomfort.

### FEATURES AND BENEFITS

- Automatic needle lock-out; appearance of red stripes confirms that the needle has been used and is ready to be disposed of safely
- Reduces user discomfort due to thinner needle (31G)
- Viewing window allows easy visibility to needle tip during priming to prevent unnecessary drug waste
- Wing-shaped design for comfortable and secure grip
- Not made with natural rubber latex



MFR#	DESCRIPTION	SPECS	QTY
16-N5MMPA	Safety Pen Needles, Viewing Window, Sterile	31G x 3/16 in	100 ea/bx, 20 bx/cs
16-N6MMPA	Safety Pen Needles, Viewing Window, Sterile	31G x 1/4 in	100 ea/bx, 20 bx/cs
16-N8MMPA	Safety Pen Needles, Viewing Window, Sterile	31G x 5/16 in	100 ea/bx, 20 bx/cs

### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.