


## LiquiBand® XL skin closure systems

McKesson LiquiBand® XL is an acrylate adhesive and mesh product that is used to hold closed easily approximated skin edges of wounds from surgical incisions up to 20 cm in length. Acts as an effective microbial barrier as long as the film remains intact.

### Features and benefits

- Provides safe, secure, and effective wound closure
- Accelerant inside applicator provides consistent dry time of 60 seconds or less
- Provides durable and effective wound closure for up to 14 days with proper after care
- Device is water-resistant once polymerized
- Simple, one-handed priming with signature winged applicator
- Wide foam tip for easy, effective adhesive application
- For use on surgical incisions, including incisions from minimally invasive surgery, and simple, thoroughly cleansed trauma induced lacerations
- Not made with natural rubber latex



### ORDERING INFORMATION

MFR #	DESCRIPTION	CONTENTS	QTY
122-LBXL	LiquiBand XL Skin Closure Systems	2 Skin Closure Systems, which each contain: • One self-adhering 8.66 in x 1.57 in mesh strip patch • One 2.7 mL 2-octyl cyanoacrylate topical skin adhesive applicator	1 system/pk, 2 pks/bx, 6 bx/cs

#### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

#### Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.