

Steam Sterilization Indicator Tape

McKesson steam sterilization indicator tape is used in securing wrapped instrument sets before, during and after sterilization. Each roll is individually packaged to ensure quality at the time of use. The blue tape features a more aggressive adhesive and is for disposable wraps only. Choose McKesson for lead-free indicator tape.

Features and Benefits

- For use in securing wrapped instrument sets before, during, and after sterilization
- Each roll of McKesson Brand tape is individually packaged to ensure quality at time of use
- Lead-free ink means no hazardous waste issues
- Color change indicator verifies exposure to sterilization process
- For use with disposable wraps only
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	SIZE	QTY
73-BT36	Sterilization Indicator Tape, Steam, Blue	1 in x 60 yds	18 ea/cs
73-ST36	Sterilization Indicator Tape, Steam	1 in x 60 yds	18 ea/cs
73-ST48	Sterilization Indicator Tape, Steam	3/4 in x 60 yds	18 ea/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908