

## Multi-Enzymatic Cleanser, Foaming

McKesson Multi-Enzymatic Cleanser Foaming Spray is specially formulated to safely remove blood, fat, tissue and body fluids from surgical instruments and equipment. The unique chemical properties penetrate and solubilize blood and fats. The wetting agents will keep instruments moist for hours, helping promote safer and easier processing in decontamination. It helps minimize handling, thus reducing exposure to potentially infectious organic contaminants. The pre-spray aids in the thorough and complete cleaning of equipment and instruments.

### Features and Benefits

- Specially designed to cling to heavy soil and prevent bioburden from drying on surgical instruments
- Non-corrosive formula helps soften and dissolve blood, fat, tissue, and other body fluids safely and effectively
- Wetting agents keep instruments moist for hours, minimizing handling and reducing exposure to potentially infectious organic contaminants
- Neutral pH formulation is safe for glass, plastic, rubber and metals as well as other surgical instruments
- When used properly, helps eliminate manual washing, instrument staining and damage
- Pleasant scent
- Phosphate-free, Chlorine-free, biodegradable



### ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Multi-Enzymatic Cleanser Foaming, 1 qt	53-28517	12/cs

#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.