

TRANZPORT HOOD

IMPORTANT PRODUCT INFORMATION—PLEASE READ BEFORE USING



INTENDED USE: To reduce the risk of the wearer transmitting fluids (saliva and mucous) from the facial area by spitting, sneezing or coughing. Improper use can result in serious injury or death. **WARNING:** Improper use may cause asphyxiation, suffocation or drowning in one's own fluids.



INSTRUCTIONS: Open bag and remove Hood. Place hood over the head of the person so the mesh fabric is covering the eye area (top) of the head, and the filtration fabric is covering the lower half of the head (below nostrils).

For the best fit, the center elastic will go under the nose and over the ears. For better protection, the elastic may be placed above the nostrils.

Carefully push the white plastic Secure-Lock tab down toward the top of the head while holding the top of the mesh fabric. This should take the slack out of the top and help secure the Hood in position. **DO NOT** push so tightly as to be uncomfortable or impair vision of the wearer.

CONDITIONS FOR USE:

- Wearer must be under control and restrained.
- Wearer must be under constant visual supervision and should **NOT** be left unattended.
- **DO NOT USE** on anyone that is vomiting, having difficulty breathing, or is bleeding profusely from the mouth or nose area.
- Remove all jewelry and eyewear before application.
- The Hood is one size fits all. If there is difficulty applying due to a large size head, discontinue use.

INVENTORY ID:

TranZport Hood, each NIK-7000
Dispenser Box of 40 NIK-7000 Box

WEIGHT: Each 1 oz.
 Box 3 lbs.

DESCRIPTION: One size fits all.

SAFETY CHECK: Including, but not limited to, fabric, mesh, elastic and adjustable tab.

CLEANING AND MAINTAINING: Dispose after use.

STORAGE: Store in a dry area to prevent mildew.

Rev 2/20/05