

Drive DeVilbiss Conserver “Tool Box”

UNIT	Specifications	Features
 Bonsai® Velocity OM-812	- Delivers up to 6:1 saving ratio up to 40 breaths per minute 1-6 uniform pulse settings, 2 CF settings preset @ 2 & 4 LPM - Operating pressure range: 500-3000 PSI 2 year warranty	- Pneumatic, requires no batteries - Quiet, compact and lightweight - Easy-to-operate single selector knob - Oxygen delivery in the first half of the inspiratory cycle
 Evolution™ OM-900	- Delivers a 5:1 savings ratio up to 40 breaths per minute - Operating pressure range: 500-3000 PSI - Weighs only 14.9 ounces with batteries - Minimum 2 year battery - utilizing 2 AA batteries 1-7 uniform pulse settings, CF preset at 2 LPM	- Electronic visual light with every breath - Auto shut-off feature conserves battery life - Larger bolus sizes, increased trigger sensitivity - Oxygen delivery in the first half of the inspiratory cycle
 Evolution™ Motion OM-900M	- Delivers a 5:1 savings ratio up to 40 breaths per minute - Operating pressure range: 500-3000 PSI - Weighs only 14.9 ounces with batteries 2 year battery life* - utilizing 2 AA batteries 1-6 uniform pulse settings, CF preset at 2 LPM	- Electronic, visual light, increased trigger sensitivity - Rest and Active (2) setting programmable - Senses patient motion and automatically adjusts the setting to match oxygen delivery to activity - Delivers oxygen in the first half of the inspiratory cycle
 SmartDose® CTOX-MN02	- Delivers up to 5:1 savings ratio up to 40 breaths per minute - Operating pressure range: 500 – 3000 PSI - Weighs just 15.5 ounces with batteries - Minimum 1 year battery life* - utilizing 2 AA batteries - Auto shut-off feature conserves battery life 1-5 uniform pulse settings, CF preset at 2 LPM - Increases 2 settings in “Active Mode”	- Electronic, Auto-adjusting, visual light increased trigger sensitivity - Continuously monitors breathing patterns and adjusts oxygen dose to match activity level - Allows delivered dose to auto-adjust up 1-2 settings during exertion to a maximum pulse dose setting of 7 - Achieves maximum conservation efficiency by automatically adjusting the dose back down when less oxygen is needed - Delivers oxygen in the first 2/3 of the inspiratory cycle
 Oxymizer® P-224 Oxymizer® O-224	- Continuous flow, disposable conserver - Provides equivalent doses of oxygen using lower liter flows up to 15 LPM of continuous flow - Delivers up to a 4:1 savings ratio - May eliminate the need to upgrade from a 5 LPM concentrator to a 10 LPM	- Provides a comfortable alternative to a mask, allowing patients to eat, drink and talk - Compatible with compressed gas cylinders, concentrators and liquid oxygen - Proven efficacy for use at rest, sleep and during exercise - Pendant model reservoir can be worn underneath clothing - Facilitates the delivery of continuous high-flow oxygen in a homecare, hospital, hospice or long-term care setting