

Medline Temperature Cable Cheat Sheet



Find the right connection.

To determine the appropriate cable, compare the connection point where the Foley catheter's temperature sensing cable attaches to the monitor with the photos listed below.

DYND110101

¼" (0.6 cm) Adapter, 12 ft/
GE Dual Temp Probes



DYND110102

2-Prong Adapter, 12 ft/
Philips



DYND110103

7-Prong Adapter, 12 ft/
Siemens



Questions? Contact your Medline Representative or visit medline.com

Medline Industries, LP | Three Lakes Drive, Northfield, IL 60093 | 1-800-MEDLINE (633-5463)

©2025 Medline Industries, LP. All rights reserved. Medline is a registered trademark of Medline Industries, LP. MKT19W7586750 / e24952 / 89