

## Multi-Use Accessory Bag

HCPCS:



- Designed with comfort and customization in mind, wear the bag over your shoulder or secure it to your mobility device using the easy-to-adjust fastening clip. The front pocket has a magnetic button closure to access items like keys and electronic devices, and a large interior pocket with O-ring zipper closure. Versatile and functional, the bag is constructed of lightweight nylon material and a reinforced back foundation for added durability.
- Back strap attachment with secure, adjustable clip
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Front pocket with magnetic closure for added security of belongings
- Reinforced back foundation maintains bag shape for added support
- Interior small pockets and zipper pocket
- Fully-lined pockets and interior
- Removable, optional cross-body carrying strap
- Treated with water resistant technology
- Spot clean with wet cloth

| PRODUCTS              | Description                         |
|-----------------------|-------------------------------------|
| RTL10255NV            | Univ Mobility Bag, Small, Navy, 1EA |
| RTL10255PK            | Univ Mobility Bag, Small, Pink, 1EA |
| SPECIFICATIONS        | Description                         |
| Weight Capacity       | 3 lb                                |
| Dimensions            | 10" (L) x 6.5" (W) x 3" (D)         |
| Retail Box Dimensions | 11" (L) x 6.5" (W) x 1.4" (D)       |
| Weight                | 0.7 lb                              |