

DOPPLERS

McKesson LUMEON® SERIES Dopplers enable medical professionals to perform accurate diagnostic assessments for a number of common procedures. The LCD display provides easy viewing of results, while the ultrasound transducers help to locate fetal and vascular targets with ergonomic comfort and excellent maneuverability. Advanced vascular system includes blood pressure cuffs to monitor a patient's cardiovascular health efficiently and comfortably. Our Lumeon Dopplers work well with probes that range from 2 to 8 MHz.

FEATURE AND BENEFITS

- Interchangeable probes to adapt to different diagnostic procedures
- Signal strength and diagnostic indicators for optimal performance
- Extended battery life and one-year warranty available
- Easy access to controls for quick viewing and reading
- Not made with natural rubber latex



MFR# 1153

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.



MFR# 1154



MFR# 1150



MFR# 1162

DESCRIPTION

DESCRIPTION	MFR#	QTY
Doppler ABI/Advanced Vascular System	1150	1ea/cs
Doppler ABI Kit	1153	1ea/cs
Doppler Handset with Display/without Probe	1154	1ea/cs
Doppler Handset without Display/without Probe	1155	1ea/cs
Doppler Handset with Display and 2 MHz Early Fetal Heart Obstetrical Probe	1156	1ea/cs
Doppler Handset with Display and 3 MHz Early Fetal Heart Obstetrical Probe	1157	1ea/cs
Doppler Handset without Display and 4 MHz Wide Beam Vascular Heart Obstetrical Probe	1158	1ea/cs
Doppler Handset without Display and 5 MHz Moderately Deep Vascular Probe	1159	1ea/cs
Doppler Handset without Display and 8 MHz Peripheral/Superficial Vascular Probe	1160	1ea/cs
Obstetrical Probe, 2 MHz	1162	1ea/cs
Obstetrical Probe, 3 MHz	1163	1ea/cs
Vascular Probe, 4 MHz	1164	1ea/cs
Vascular Probe, 5 MHz	1165	1ea/cs
Vascular Probe, 8 MHz	1166	1ea/cs