

Rad-G™ Pulse Oximeter


A rugged handheld device powered by clinically proven Masimo SET® pulse oximetry



Strong external casing built to withstand drops from 6 ft. onto a hard surface



Light weight and slim profile allow for easy transport and storage



Internal rechargeable battery provides 24 hours of operational use



Reliable Monitoring

Hospital-grade Masimo SET® Measure-through Motion and Low Perfusion™ pulse oximetry*

Informative Display

Display of pleth waveform for additional information to use in patient assessment

Alarm Management

Configurable alarms to offer insight into patient status changes that may require intervention



Available with Masimo single-patient-use and reusable sensors, Rad-G monitors:



Clinically Proven Technology Used in Leading Hospitals Worldwide

per year around the world¹

hospitals according to the 2020-2021
U.S. News & World Report
Best Hospitals Honor Roll²

have shown that Masimo SET®
outperforms other pulse
oximetry technologies³

ACCURACY ARMS⁴

ACCURACY ARMS ⁴		ENVIRONMENTAL	
Oxygen Saturation (SpO2)	70–100%	Operating Temperature	While battery is charging: ⁵ 0–40° C (32–104° F) While battery is NOT charging: ⁶ 0–50° C (32–122° F)
No Motion Adults/Pediatrics/Infants	2%		
Motion Adults/Pediatrics/Infants	3%		
Low Perfusion Adults/Pediatrics/Infants	2%		
Pulse Rate (PR)	25–240 bpm	BATTERY	
No Motion Adults/Pediatrics/Infants	3 bpm	Battery Type	Lithium ion
Motion Adults/Pediatrics/Infants	5 bpm	Capacity	24 hours ⁷
Low Perfusion Adults/Pediatrics/Infants	3 bpm		
Respiration Rate (RRp)	4–70 rpm	ORDERING INFORMATION	
No Motion Adults/Pediatrics (>2 years of age)	3 rpm ARMS ⁴ , ± 1 rpm mean error	Rad-G Pulse Oximeter Kit with Sensor	PN 9847
		Rad-G Pulse Oximeter Kit	PN 9849
		Rad-G Reusable Sensor	PN 4325
		Rad-G Reusable YI Sensor	PN 4653
		RD SET® G15-05 Patient Cable	PN 4773
		RD SET G15-12 Patient Cable	PN 4774
PHYSICAL CHARACTERISTICS			
Weight	0.27 kg. (0.59 lbs.)		
Dimensions	7.4 cm x 19.8 cm x 2.5 cm (2.9" x 7.8" x 1.0")		

Not licensed for sale in Canada.

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

